


DISSECTING KITS




MANUFACTURER & EXPORTER

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-001

Dissecting Set of 9 Pcs.

- Dissecting Forceps 11.5cm Sharp Point.
- Teasing needle Sharp STR.
- Teasing Needle Sharp Angular.
- B.P Handle No.4.
- Blade No. 22.
- Thumb Forcep T.O.E 13cm.
- Dissecting Scissor STR 11.5cm.
- Dissecting Scissor STR 13cm.
- Zipper Case Inside Red Velvet.



TC17-002

Dissecting Set of 10 Pcs.

- Grooved Director with Tongue Tie 14cm.
- Probe Double 14cm.
- Mayo Hegar Needle Holder 14cm.
- B.P Handle No.3
- B.P Handle No.4
- Dressing Scissor 14cm S/B STR.
- Dressing Scissor 14cm S/B CVD.
- Dressing Forcep 1x2 teeth 14cm.
- Dressing Forcep Serrated 14cm.
- Zipper Case Inside Red Velvet.



TC17-003

Dissecting Set of 10 Pcs.

- Probe Double 13cm.
- 2PCS. Scalpel S.S.
- Grooved Director With Tongue Tie 14cm
- Mayo Hegar Needle Holder 14cm.
- Dissecting Forceps 14cm 1x2 Teeth.
- Dressing Forceps 14cm Serrated.
- Dressing Scissor 14cm S/B STR.
- Dressing Scissor 14cm S/B CVD.
- Zipper Case Inside Red Velvet.



TC17-004

Dissecting Set of 10 Pcs.

- Blow Pipe.
- 3PCS Scalpel.
- Dissecting Forceps 14cm.
- Grooved Director With Tongue Tie 14cm
- Triple Chain Hook.
- Dressing Scissor Probe Pointed 14cm.
- Probe Double 14 cm.
- Zipper Case Inside Red Velvet.



Dissecting Set


TC17-005 Containing 9 Instruments

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico


Combinazioni per diagnostica



17-006

Basic Dissecting Set

A Collection of Instruments For Introductory Dissection vinyl pouch
Scalpel 6" Scissors 4 1/2" Forceps student grade 4 1/2"
Straight needle in Wooden Handle Bent needle in wooden handle
Medicine Droppers.



17-007

Dissecting Set

Beginners For Elementary Junior High Or High School Life Science
Vinyle pouch Scissors Forceps student /Scalpel with Replaceable
blade Straight needle in wooden handle Bent needle in wooden
handle Medicine dropper



17-008

General Biology Dissecting Set


Designed for work in general biology good Quality instruments
scissor student Forceps fine point 4 1/2" S/S
Scalpel 1 1/2" Blade (with G pin S/S)
Straight Needle with wooden handle Medicine Dropper
single fold leatherette case.



17-009

Fine Dissecting Set


Designed for Freshman Courses in Biology
Scissors Fine point 4 1/2" Forceps fine point with g pin 4 1/2"
Scalpel replaceable Blade Teasing Needle Str
Teasing Needle Cvd Medicine Dropper
Single Fold Leatherette Case.



17-010

Intermediate Dissecting Set


Contains high quality instruments recommended for advance
biology courses . Scissor surgical 4 1/2" Forceps fine point 4 1/2 "
Surgeon scalpels (with # 22 blade) Teasing needle plastic handle Str
Teasing needle Cvd Plastic handle Mall probe and seeker
Medicine dropper Single fold leatherette case



17-011

Tissue Culture Dissecting Set

Features Stainless Steel Instrument in leatherette case
Scissors M/Dissecting Iris 4 1/2 " Forceps watchmaker 4 1/4" Str
Forceps watchmake 4 1/4" Cvd Scalpel micro dissecting 5"
18x3 mm blade Teasing Needle Str S/S Teasing Needle Cvd S/S
Single Fold leatherette case



TC17-012

Anatomy Dissecting set

A complete set Detailed Dissection for more advanced work in zoology and comparative anatomy .
Surgical scissors 5 1/2". Fine point scissors 4 1/2" .
Forceps fine point Cvd 4 1/2" Thumb forceps superior 5" .
Scalpel 1 1/2" Blade .
Cartilage knife/Flexible Probe. Triple chain Hook/Mall probe.
Needle wooden handle Str. Needle wooden handle Str
Medicine Dropper. Double fold leatherette case



TC17-013

Graduate Dissecting Set

Design for the biology graduate students
Surgical scissors 5 1/2". Fine point scissors 4 1/2" .
Forceps fine point Str 4 1/2" Forceps fine point Cvd 4 1/2" .
Tripple chain hook. cartilage knife scalpel stain finish 1 1/2" blade
forceps cartilage 5" 2 needle holder with straight needle .
Flexible probes. Probes with hook. Probe & seeker
Dropping Pipe Double Fold Leatherette case



TC17-014

Micro Dissecting Set

Feature stainless steel Instruments in Leatherette case.
Scissors M/Dissecting Iris 4 1/2" Scissors sharp point Iris 4 1/2" .
Forceps Micro Dissecting 4 1/4" Straight Forceps micro dissecting 4" .
Curved Scalpel micro dissecting 18 x 3 mm blade
scalpel micro dissecting 30 x 3 mm blade Dissecting needle Str .
Dissecting needle Cvd. Dissecting needle Full-Spear point
Dissecting needle Half-Spear point In double folding case.



TC17-015

Large Animal set

This set provide all the Necessary instruments Dissection of large Animal such as Cat, Shark and Rabbits Scissor surgical 5 1/2" .
D/Scissors superior 4 1/2" Scissors fine point Str 4 1/2" .
Forceps large Tissue 1x2 teeth 6" Forceps utility Cvd point 6" .
Forceps fine point str 4 1/2" Scalpel screw Lock with # 22 blade .
Knife cartilage 6 1/2" . Probe sharp point angular. Section lifter S/S
Teasing needle Str Wooden handle. Teasing Needle Cvd Wooden handle
Medicine dropper.

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico


Combinazioni per diagnostica



TC17-016

New Instructor's Dissecting set

Contain quality stainless steel Instrument Enclosed in an attractive zippered case Scissor surgical 4 1/2" scissors fine point 4 1/2" Forceps fine point Str. Forcep fine point Cvd. Mosquito forceps Cvd 5" box lock. One piece scalpel 1 1/2" blade Surgeon handle # 4. Mall probe and seeker 2 plastic handle str needle also 2 blade # 22 for each set.



TC17-017

Dissecting Set- Advanced

Stainless steel, consisting of 12 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs. Artery Forceps - 1 pc. Tissue Forceps - 1 pc. Seeker - 1 pc. Scalpel Handle w/ Blade - 1 pc. Bone Cutter - 1 pc. Dental Probe & Seeker - 1 pc.



TC17-018

Dissecting Kit

Kit includes:

1 4 1/2" Dissecting forceps with curved and serrated tips
1 Curved dissecting needle 1 Glass pipette with rubber bulb
1 4 1/2" Dissecting forceps straight and serrated tips
1 Dissecting scalpel, 6" long, 1 1/2" cutting edge.
1 4 1/2" Dissecting scissors, sharp, straight points
1 Plastic ruler 6", 1/16" div., millimeter and centimeter



TC17-019

Dissecting Set - Pro

Stainless steel, consisting of 20 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc. B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Ruler 6" - 1 pc. Magnifying Glass - 1 pc. Glass Slide - 1 pc. Chain and Hook - 1 pc. Cover Slips - 1 pc. Blow Pipe - 1 pc. Needle Wooden Handle - 1 pc. Dropper Glass - 1 pc. Brush - 1 pc.



TC17-020

Dissecting Set - Popular

Stainless steel, consisting of 14 instruments for accurate dissection. Complete in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc. B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Magnifying Glass - 1 pc. Chain and Hook - 1 pc. Blow Pipe - 1 pc.



TC17-021

Dissecting Set - Senior

Stainless steel, consisting of 9 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs. Scalpel - 1 pc. Scalpel Handle w/ Blade - 1 pc. Ruler 6" - 1 pc.

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-022

Dissecting Set - Student

Stainless steel, consisting of 9 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

- Scissors - 2 pcs. Forceps - 2 pcs.
- Needle Wooden Handle- 2 pcs.
- Scalpel - 1 pc. Scalpel Handle w/ Blade - 1 pc.
- Dropper Glass - 1 pc.



TC17-023

Dissecting Set - Beginner

Stainless steel, consisting of 7 instruments for accurate dissection. Complete in zipper leather case.

Kit consists of: Scissors - 2 pcs.

Forceps - 2 pcs. Seeker - 1 pc.

B.P. Handle w/ Blade - 1 pc.

Needle Wooden Handle - 1pc



TC17-024

MICROSCOPIC DISSECTING SET



TC17-025

Dissecting-kits

11 Pcs



17-026

Biology Dissecting Tools Set



TC17-027

Anatomical Dissecting Set


13 parts filled. in leatherette zipper pocket
Content

- 1 straight dissecting needle with plastic handle, 140 mm
- 1 bent dissecting needle with plastic handle, 140 mm
- 1 lancet with plastic handle, 140 mm
- 1 bandage scissors, pointed/blunt, 130 mm
- 1 straight microscope scissors, 115 mm
- 1 bent microscope scissors, 105 mm
- 1 straight pointed tweezers, 105 mm
- 1 straight blunt tweezers, 115 mm
- 1 pointed scalpel, 130 mm, cutting edge 17 mm
- 1 convex scalpel, 130 mm, cutting edge 17 mm
- 1 pointed scalpel, 155 mm, cutting edge 40 mm
- 1 section lifter, 160 mm, width 10 mm
- 1 directional probe, 140 mm



TC17-028

DE-SET1



Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico


Combinazioni per diagnostica



TC17-029

Dissecting Set - Economy

Stainless steel, consisting of 7 instruments for accurate dissection. Supplied in leather storage case. Kit consists of: Scissors - 1 pc. Forceps - 1 pc. Needles - 2 pcs. B.P. Handle w/ Blade - 1 pc. Dropper Glass - 1pc. Ruler 6" - 1 pc.



TC17-030

Basic Dissecting Kit



TC17-031

Basic Dissecting Tools Set



TC17-032

Deluxe Instructor's Dissecting Set



TC17-033

Podiatry Kit



TC17-034


Anatomy Dissecting Set



TC17-035

Micro Dissecting Set

Features stainless-steel instruments in Double Leatherette Case



TC17-036
Standard Dissecting Set II



TC17-037

Thomas Basic Dissecting Kit

Designed for students or limited use 6 instruments
Chrome and stainless steel construction Kit includes:
6" chrome scalpel
4 1/2" Stainless steel dissecting scissors
4 1/2" Chrome-plated forceps
Straight dissecting needle Bent dissecting needle
Pipette 6" Ruler Single fold vinyl case with lining



TC17-038
General Biology Dissecting Set



TC17-039
Dissecting Set, Biology Training, Molded Plastic Case



TC17-040
Advanced Anatomy Dissecting Set



TC17-041
Advanced Anatomy Dissecting Set
with Leather Case

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-042
Teacher's Dissecting Set



TC17-043 Comparative Anatomy Dissecting Sets

Set is sold with either Double Leatherette Case or Molded Plastic Case.



TC17-044
Dissecting Set



TC17-045
MICROSCOPIC DISSECTING SET

5 parts or 7 parts filled, in wooden case Content

- 1 straight scissors, 100 mm
- 1 bent dissecting needle, 150 mm
- 1 straight dissecting needle, 150 mm
- 1 tweezers, straight point, with guide pin, 115 mm
- 1 scalpel 125 mm



TC17-046
CHEMISTS' TOOL SET
4 parts filled, in a case with press stud
Content
1 powder spatula, 150/9 mm
1 micro-powder spatula 130/4 mm
1 pointed/blunt scissors, 130 mm
1 straight pointed tweezers. 130 mm



TC17-047
Dissecting Kit

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-048
Dissecting Set



TC17-049
Dissecting Set



TC17-050
Dissecting Set



TC17-051
Dissecting Set



TC17-052
Dissecting Set



TC17-053
Dissecting Set

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-054
Dissecting Set



TC17-055
Dissecting Set




TC17-056
Dissecting Set



TC17-057
Dissecting Set



TC17-058
Dissecting Set



TC17-059
MICROSCOPIC DISSECTING SET

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-060
Deluxe Dissecting Tools Set



TC17-061
Dissecting Set



TC17-062
DISSECTING SET



TC17-063
DISSECTING SET



TC17-064
DISSECTING SET



TC17-065
DISSECTING SET

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-066

Advanced Biology Kit

Advanced Biology includes nine stainless steel instruments and two surgical blades tucked in a compact vinyl zippered case with red velvet lining.



TC17-067

Anatomy Dissection Kit II



TC17-068

Anatomy Dissection Kit



TC17-069


Teachers Dissection Kit

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica




TC17-074

Dissection Kit w/ Screw Lock Blade



TC17-075

Precision Dissection kit
Precision Dissection kit - plastic case



TC17-076

Dissection Kit - Basic



TC17-077


Budget Dissection Kit

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico


Combinazioni per diagnostica



TC17-078
Biology Dissection Kit



TC17-079
Anatomy Dissection Kit & Stethoscope Package
Med-Student first year gross anatomy dissection kit



TC17-080
Zippy Precision Kit



TC17-081
Introductory Anatomy Dissection Kit

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico


Combinazioni per diagnostica



TC17-086
Precision Dissection Kit Package with
Safety Eyewear and Gloves



TC17-087
Comprehensive Dissection Kit w/
Goggle - Kit-5



TC17-088
Entomology Dissection Kit
for insect studies.



TC17-089
Dissection Kit - Intermediate III - Kit-2

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica

	 TC17-090 Student Dissection Kit (Plastic Case)
	 TC17-091 Zippy Dissection Kit w/ Screw Lock Scalpel
	 TC17-092 Veterinary Student Dissection Kit First Year Veterinary Student Anatomy Training Kit
	 TC17-093 Podiatry Students Training Kit Podiatrist Student Dissecting Kit

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-094

Zippy Dissection Kit Package With Nitril Gloves and Safety Eyewear

Ideal dissection package includes Zippy dissection set, eyewear, and pack of nitrile gloves




TC17-095

Comprehensive Economy Dissection Kit

Comprehensive Dissecting kit includes 13 instruments and 2 packs of blades to set any student up for a science lab.

All tools come neatly packed away in a snap button case.



TC17-096

Student Anatomy Dissection Kit - Economy



TC17-097

Invertebrate Dissection Kit

Includes 8 tools and 2 blades in a sleek vinyl zippered case.

All instruments are stainless steel unless otherwise noted.

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica

	 TC17-098 Budget Entomology Kit Budget Entomology dissection kit includes the basic tools needed for a classroom dissection.
	 TC17-099 General Biology Dissection. Kit Basic General Biology Dissecting Kit
	 TC17-100 Med Student Anatomy Dissection Kit Anatomy dissection kit for first year medical students!
	 TC17-101 Dissection Kit Standard

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-102
Med Student Anatomy Dissection Kit
Anatomy Dissection Kit for
first year medical students



TC17-103
Dissection Kit Deluxe
Deluxe Dissecting Kit In A
Zippy Case



TC17-104
Botany Dissection Kit - Economy
Budget Botany Dissection kit



TC17-105
DR Botany Dissection Kit
Botany Dissection Kit

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-106
Dissecting Set



TC17-107
Dissecting Set



TC17-108
Dissecting Set



TC17-109
Dissecting Set



TC17-110
Student Dissecting Set I



TC17-111
Intermediate Dissecting Set I

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica



TC17-112



TC17-113



TC17-114



TC17-115



TC17-116
Dissection Premium Kit



TC17-117
Dissection Premium Kit



TC17-118
Set of 6 Pcs

Double Spoons, Length 150 mm, L x I spoon 1: 11/22 mm, L x I spoon 2: 20/25 mm
Double Spoons, Length 180 mm, L x I spoon 1: 20/27 mm, L x I spoon 2: 24/32 mm
Double Spoons, Length 210 mm, L x I spoon 1: 20/28 mm, L x I spoon 2: 27/35 mm
Double Spoons, Length 250 mm, L x I spoon 1: 24/32 mm, L x I spoon 2: 30/40 mm
Double Spoons, Length 300 mm, L x I spoon 1: 29/40 mm, L x I spoon 2: 35/45 mm
Double Spoons, Length 160 mm, L x I spoon 1: 22/35 mm, L x I spoon 2: 30/45 mm



TC17-119
Set of 2 Pcs

Double spoons, length 160 mm, L x I spoon 1: 22/35 mm, L x I spoon 2: 30/45 mm
Double spoons, length 230 mm, L x I spoon 1: 30/60 mm, L x I spoon 2: 38/75 mm



TC17-120
Set of 8 Pcs



Double spatulas, spatula width 15 mm, spatula length 70 mm, length 300 mm
Double microspatula, straight edges, spatula width 6 mm, spatula length 40 mm, length 130 mm
Double microspatula, straight edges, spatula width 6 mm, spatula length 40 mm, length 150 mm
Double microspatula, straight edges, spatula width 6 mm, spatula length 50 mm, length 185 mm
Double microspatula, straight edges, spatula width 6 mm, spatula length 55 mm, length 210 mm
Double microspatula, trapezoid shape, spatula width 5 mm, spatula length 28 mm, length 130 mm
Double microspatula, trapezoid shape, spatula width 5 mm, spatula length 30 mm, length 185 mm
Double microspatula, trapezoid shape, spatula width 5 mm, spatula length 30 mm, length 210 mm




TC17-121
Micro-spatula set
with plastic handle, stainless steel
Set of 4 parts
Content:
each 1 spatula 3, 4, 5, 6mm wide
Length in mm 150



TC17-122
Micro-powder Spatula Set
with plastic handle, stainless steel
Set of 4 parts
Content:
each 1 spatula 3, 4, 5, 6mm wide
Length in mm 150



TC17-123
Set of 2 Pcs
Excavator with button
From Hammacher. Stainless steel Remanit 4024. Autoclavable.
Shovel with button. Length 190 mm. Excavator dimensions: 35/70
Shovel with button. Length 205 mm. Excavator dimensions: 30/63



TC17-124
Set of 6 Pcs
Micro Double Spatula stainless steel
flexible

length mm	width of blade mm
100	3 2 g
130	3 2 g
150	3 4 g
185	3 4 g
150	6 5 g
185	6 6 g



TC17-125
Set of 6 Pcs
Double Spatula stainless steel
flexible

length mm	width of blade mm
130	9 6 g
150	9 7 g
185	9 8 g
210	11 15 g
250	11 17 g
300	16 42 g



TC17-126
Set of 3 Pcs
Double Spatula for Mortars
stainless steel
width of blade 20 mm rigid

length mm	weight g
150	29 g
210	38 g
250	49 g



TC17-127
Set of 3 Pcs
Powder Spatula stainless steel
width of blade 9 mm

length mm	weight g
150	8 g
185	8 g
210	10 g



TC17-128
Set of 4 Pcs
Double Spatula bent
stainless steel curved

length mm	width of blade mm
150	9 7 g
185	10 9 g
210	12 16 g
250	12 18 g



TC17-129
Set of 3 Pcs
Micro Spatula Spoon
stainless steel rigid

length mm	weight g
150	6 g
185	7 g
210	8 g



TC17-130

Set of 4 Pcs

Spatula Spoon stainless steel
bent spoon width of spoon 19 mm

length mm

150	19 g
185	23 g
210	28 g
250	33 g



TC17-131

Set of 4 Pcs

Spatula Spoon stainless steel
rigid broad

length mm

length mm	width of blade mm	
120	23	15 g
150	28	23 g
180	30	35 g
210	30	42 g



TC17-132

Set of 3 Pcs

Spatula Spoon stainless steel
width of spoon 29 mm

length mm

150	20 g
180	30 g
210	37 g



TC17-133

Set of 6 Pcs

Double Spoon stainless steel
rigid round handle

length mm

length mm	width of spoon mm	
120	17 / 25	16 g
150	17 / 25	19 g
180	17 / 25	29 g
210	22 / 29	34 g
250	29 / 35	48 g
300	29 / 35	61 g



TC17-134

Chemist's Spoon
stainless steel
length 150 mm
width 16 mm
12 g



TC17-135

Spoon, open
stainless steel
length 200 mm
width 25 mm
29 g



17-136
Set of 2 Pcs
Sodium Spoon
copper basket
basket Ø 25 mm
length 180 mm
Single 16 g
double 26 g



TC17-137
Deflagrating Spoon
cup stainless steel
cup Ø 16 mm
length 450 mm
14 g



TC17-138
Ladle
stainless steel
cup Ø 75 mm
length 275 mm
73 g



TC17-139
Set of 7 Pcs
Spatula stainless steel
highly flexible polished, wooden handle
length of blade mm length mm
75 155 16 g
100 190 19 g
130 230 29 g
150 255 34 g
200 320 48 g
250 360 61 g
300 405 61 g



TC17-140
Set of 5 Pcs
Spatula stainless steel highly
flexible plastic handle
autoclavable
length of blade mm length mm
75 160 19 g
100 200 26 g
130 240 31 g
150 260 35 g
200 330 63 g



TC17-141
Spatula stainless steel
blade flexible
length of blade 50 mm
width of blade mm
15 23 g



Wazirbad Road, Ugoki, Sialkot-51310 PAKISTAN

Phone: +92-52-3561347-3553588

Phone: +92-52-3561202

Cell: +92-300-6139615

Cell: +92-322-6139615

Email: info@tcmedicline.com

Email: exporttcmedicline@gmail.com