

DIAGNOSTIC PERIODONTAL & RESTORATIVE INSTRUMENTS

DENTAL HAND INSTRUMENTS

The types of hand instruments can be categorized based on their function:

- Examining the mouth and teeth
- Scaling or professional cleaning of teeth
- Cutting teeth and removing cavities
- Placing and condensing filling materials
- Carving and finishing dental fillings
- Miscellaneous types
- Instruments used for examining the mouth and teeth

Mouth mirrors are used to examine the condition of the teeth and mouth, as well as to retract the tongue and cheeks upon examination.

Probes— There are many probes available but the main types are the straight probe, Briault probe and periodontal probe.

The straight probe is sharp and it is used to examine cavities in the dentine region as well as to check the margins of restorations or fillings.

The Briault probe is also sharp and it is used to detect the junction between the enamel and dentine. Briault can also be used to detect tartar hidden in the gum pockets.

A periodontal probe can be used to measure the depth of the gum pockets and for other measurements. Instruments for cutting and scraping tissues and deposits

Scalers remove tartar (calculus) on the tooth surface as well as those hidden in the gum pockets. They can also be used to remove other surface deposits and temporary crowns.

Excavators are used to remove softened dentine and temporary fillings. The back of the blade can also be used for placing linings between filling materials. Sometimes excavators can be used for carving metal or amalgam fillings.


Chisels are used to remove unsupported enamel during cavity preparation. Unsupported enamel can be harmful for it can compromise the longevity of the filling if not removed.

Hatchets and hoes have similar purposes to chisels. Instruments used for placing and condensing filling materials. Plastic instruments convey and shape filling materials and are usually flat blades that do not use heavy pressure.

Some instruments have rounded ends for shaping and burnishing the surfaces of fillings. Plastic instruments are usually made of stainless steel for most materials though there are also thin Teflon-coated or titanium nitride instruments which are non-stick for composite fillings. **Condensers or pluggers** are used for compressing and forming filling materials. Heavy pressure is applied on use therefore these instruments are usually used for metal or amalgam fillings. Instruments used for carving and finishing fillings. Carving and finishing instruments have sharp or semi-sharp blades of various shapes. They are used to carve materials by cutting and scraping.



Symptoms of Gum and Dental Problems



MANUFACTURER & EXPORTER

DIAGNOSTIC PERIODONTAL & RESTORATIVE INSTRUMENTS

VALUED BUYERS.

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours. We regularly manufacturing below articles in showing handles & pattern A to L for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in materials, sizes, shapes and patterns.



A-Hollow Handle Ø 10 mm



B-Hollow Handle Ø 10 mm



C-Hollow Handle Ø 10 mm



D-Hollow Handle Ø 10 mm



E-Hollow Handle Ø 10 mm



F-Hollow Handle Ø 10 mm



G-Hollow Handle Ø 10 mm



H-Hollow Handle Ø 6.4 mm



I-Hollow Handle Ø 8 mm



J-Hollow Handle Ø 10 mm



K-Soild Handle Ø 6 mm



L-Soild Handle Ø 6 mm

MANUFACTURER & EXPORTER

зеркал
Espejos de la boca
Mundspiegel

дзеркал
Miroirs à la bouche
مرآة الفم

Plain



TC27-001 Fig. 3 Ø20mm
TC27-002 Fig. 4 Ø22mm
TC27-003 Fig. 5 Ø24mm
TC27-004 Fig. 6 Ø26mm

Cone-Socket Plain



TC27-005 Fig. 3 Ø20mm
TC27-006 Fig. 4 Ø22mm
TC27-007 Fig. 5 Ø24mm

Magnifying



TC27-008 Fig. 3 Ø20 mm
TC27-009 Fig. 4 Ø22 mm
TC27-010 Fig. 5 Ø24 mm

Cone-Socket Magnifying




TC27-011 Fig. 3 Ø20mm
TC27-012 Fig. 4 Ø22mm
TC27-013 Fig. 5 Ø24mm

Rhodium



TC27-014 Fig. 3 Ø20mm
TC27-015 Fig. 4 Ø22mm
TC27-016 Fig. 5 Ø24 mm

Cone-Socket Rhodium



TC27-017 Fig. 3 Ø20mm
TC27-018 Fig. 4 Ø22mm
TC27-019 Fig. 5 Ø24mm

Front Surface



TC27-020 Fig. 3 Ø20mm
TC27-021 Fig. 4 Ø22mm
TC27-022 Fig. 5 Ø24mm
TC27-023 Fig. 6 Ø26mm

Cone-Socket Front Surface



TC27-024 Fig. 3 Ø20mm
TC27-025 Fig. 4 Ø22mm
TC27-026 Fig. 5 Ø24mm

Ручки для зеркал

Manches pour miroirs à bouche

Mundspiegelgriffe

Ручки для дзеркал

Mangos para espejos de boca

مقابض لمرآة الفم



TC 27-027
Hollow Mirror
Handle Round
Ø 10 mm



TC 27-028
Hollow Mirror
Hexagonal
Handle
Ø 8 mm



TC 27-029
Hollow
Ø 8 mm



TC 27-030
Hollow Mirror
Hexagonal
Handle
Ø 10 mm



TC 27-031
Hollow Mirror
Handle
Ø 8 mm



TC 27-032
Hollow Mirror
Handle
Ø 10 mm



TC 27-033
Hollow Mirror
Handle
Ø 8 mm



TC 27-034
Mirror Handle
Autoclavable
Material
Available in
Various Colours
Ø 9.5 mm



TC 27-035
Mirror Handle
Autoclavable
Material
Available in
Various Colours
Ø 9.5 mm



TC 27-036
Mirror Handle
Silicon Grip
Autoclavable
Available in
Various Colours



TC 27-037
Mirror Handle
Solid
Ø 6 mm



TC 27-038
Mirror Handle
Solid
Ø 6 mm



TC 27-039
Mirror Handle
Solid Octagonal
Ø 6 mm



TC 27-040
Mirror Handle
Solid
Stem with
Rule 35mm



TC 27-041
Cone-Socket
Solid Octagonal
Ø 6 mm



TC 27-042
Cone-socket
Solid
Ø 6 mm



TC 27-043
Cone-socket
Hexagonal
Handle
Ø 6 mm

Cone-Socket Mirror Handle

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقف



TC27-044
London College
15cm



TC27-045
Meriam
16cm



Semken Taylor
TC27-046 12.5cm
TC27-046A 15cm
Straight



Semken Taylor
TC27-047 12.5cm
TC27-047A 15cm
Curved



Semken Taylor
TC27-048 12.5cm
TC27-048A 15cm
1x2 Straight



Semken Taylor
TC27-049 12.5cm
TC27-049A 15cm
1x2 Curved



TC27-050
Flagg
15cm



TC27-051
Meriam
16cm



TC27-052
Amsterdam
15.5cm



TC27-053
Perry
13cm
Curved



TC27-054
Self-Locking
15cm
Angled



TC27-055
Endo-Locking
Grooved
15cm

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten


 Пінцети
 Pinzas para algodón
 ملاقط



TC 27-056 9cm
TC 27-057 12cm
TC 27-058 15cm
TC 27-059 18cm
 Taylor



TC 27-060 9cm
TC 27-061 12cm
TC 27-062 15cm
TC 27-063 18cm
 Taylor



TC 27-064 0,6mm
TC 27-065 0,9mm
 Gerald
 17.5cm



TC 27-066 0,6mm
TC 27-067 0,9mm
 Gerald
 17.5cm



TC 27-068
 Gerald
 17.5cm



TC 27-069
 Gerald
 17.5cm



TC 27-070
 15cm
 College



TC 27-071
 Meriam
 15cm



TC 27-072
 De Backy
 Atraumatic
 15cm
 Straight



TC 27-073
 Meriam Tweezers
 smooth 16cm



TC 27-074
 London College
 smooth 16cm



Adson With Teeth
TC 27-075 12cm
TC 27-076 15cm



 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



TC27-077
Perry
10.5 cm



TC 27-078
Perry
10.5 cm



TC27-079
Fig. 1
Perry
13 cm



Potts-Smith Dressing Forceps
TC27-080 18 cm
TC27-081 20 cm
TC27-082 25 cm
 straight serrated



Potts-Smith Dressing Forceps
TC27-083 18 cm
TC27-084 20 cm
TC27-085 25 cm
 curved serrated



Potts-Smith Tissue Forceps 1x2 teeth
TC27-086 18 cm
TC27-087 20 cm
TC27-088 25 cm
 straight



Potts-Smith Tissue Forceps 1x2 teeth
TC27-089 18 cm
TC27-090 20 cm
TC27-091 25 cm
 curved



De-Bakey With Tungsten Carbide (TC) 17.5cm
TC27-092



Adson
TC27-093 12cm
TC27-094 15cm



Adson 1x2
TC27-095 12cm
TC27-096 15cm



Adson Brown
TC27-097 12cm
TC27-098 15cm



TC27-099
 Crane Kaplan
 15cm



TC27-100
 Hemingway
 15cm

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пінцети
 Pinzas para algodón
 ملاقط



TC27-101
Crane Kaplan
15cm



TC27-102
Crane Kaplan
15cm



Dressing Forceps

TC27-103 10.5cm
TC27-104 11.5cm
TC27-105 13.0cm
TC27-106 14.0cm
TC27-107 16.0cm
TC27-108 18.0cm
TC27-109 20.0cm
TC27-110 25.0cm
TC27-111 30.0cm



Tissue Forceps
1x2 teeth

TC27-112 10.5cm
TC27-113 11.5cm
TC27-114 13.0cm
TC27-115 14.0cm
TC27-116 16.0cm
TC27-117 18.0cm
TC27-118 20.0cm
TC27-119 25.0cm
TC27-120 30.0cm




Tissue Forceps
2x3 Teeth

TC27-121 11.5 cm
TC27-122 13.0 cm
TC27-123 14.0 cm
TC27-124 16.0 cm
TC27-125 18.0 cm
TC27-126 20.0 cm
TC27-127 25.0 cm
TC27-128 30.0 cm



Lane Tissue Forceps




TC27-129
12.5 cm
1x2 teeth



TC27-130
15 cm
1x2 teeth



TC27-131
17.5 cm
1x2 teeth



TC27-132
15 cm
2x3 teeth



TC27-133
17.5 cm
2x3 teeth



Adson-Brown

TC27-134
12 cm
Tissue, 7:7 Teeth

TC27-135
12 cm
Tissue, 9:9 Teeth



Adson-Brown

TC27-136
12 cm
Teeth 9:9
Curved

Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

Пінцети
 Pinzas para algodón
 ملاقط

TC27-137
 15 cm
 Straight, 9:9 Teeth

TC27-138
 12 cm
 Dressing, Serrated

TC27-140
 12 cm
 dressing, serrated

TC27-141
 12 cm
 dressing,
 Micro serrated

TC27-139
 12 cm
 Tissue, 1:2 Teeth

TC27-142
 12 cm
 Tissue,
 1:2 Teeth

TC27-143
 12 cm
 Tissue, Micro
 1:2 Teeth

Brown

Ewald / Hudson

Adson /
Micro-Adson

TC27-144
 12 cm
 Dressing,
 cross serration

TC27-145
 12 cm
 cross serration-
 1:2 Teeth

TC27-146
 15 cm
 cross serration-
 1:2 Teeth

TC27-147
 Adson
 1 mm
 12 cm

TC27-148
 Adson T.C
 1,7 mm
 12 cm

TC27-149
 Adson
 1,7 mm
 15 cm

TC27-150
 Micro-Adson
 very fine
 0,8 mm
 15 cm

Adson
T.C

TC27-151
 Adson
 1x2, 1,5 mm
 12 cm

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقف



TC27-152
Micro-Adson
1x2, 1 mm
15 cm
very fine



TC 27-153
Micro-Adson
1x2, 1 mm
15 cm
curved
very fine



TC27-154
10 cm

TC 27-155
10 cm

TC 27-156
10 cm



TC27-157
10 cm
1x2 teeth

TC27-158
10 cm
1x2 teeth

TC27-159
10 cm
1x2 teeth



TC27-160
15 cm
1x2 teeth

Graefe

Graefe Tissue Forceps

Mc Indoe



TC27-161
Cushing
1,2 mm
18 cm



TC27-162 13 cm
1x2, 1,5 mm
27-163 14 cm
1x2, 1,7 mm



TC27-164
Taylor
1x2, 1,5 mm
17,5 cm



TC27-165 13 cm
1x2, 2,3 mm
27-166 14 cm
1x2, 2,3 mm



TC27-167
Micro-Cushing
1x2, 0,9 mm
17,5 cm
very fine



TC27-168
Suture Forceps
15 cm

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقف



Delicate Tweezers

TC 27-169
11.5 cm

TC27-170
13 cm

TC27-171
14.5 cm



Gerard

TC27-172
17.5 cm

TC 27-173
23 cm

TC 27-174
17.5 cm



Gerard 1x2

TC27-175
17.5 cm

TC 27-176
23 cm

TC27-177
17.5 cm



Selman
TC27-178
20 cm

Ø 4 mm



Brown

TC27-179
15 cm

TC27-180
20 cm



Strassburger Modell

TC27-181
15 cm

TC27-182
20cm

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقط



Taylor Dressing Forceps
TC 27-183 18 cm
straight serrated



Taylor Tissue Forceps
TC 27-184 18 cm
1x2 teeth, Straight



Taylor Tissue Forceps
TC 27-185 18 cm
1x2 teeth, Curved



Cushing Dressing Forceps
TC 27-186 18 cm
Serrated Straight



Cushing Dressing Forceps
TC 27-187 18 cm
Serrated Curved



Cushing Tissue Forceps
1x2 teeth
TC 27-188 20 cm
Straight



Cushing Tissue Forceps
1x2 teeth
TC 27-189 18 cm
Curved



Brophy Dressing Forceps
TC 27-190 20 cm
straight serrated



Brophy Dressing Forceps
TC 27-191 20 cm
curved serrated



Brophy Tissue Forceps
1x2 teeth
TC 27-192 20 cm
straight serrated




Brophy Tissue Forceps
1x2 teeth
TC 27-193 20 cm
curved serrated

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



Gillies Tissue Forceps
1X2 Teeth
TC27-194 15 cm
 straight serrated



Gillies Tissue Forceps
1X2 Teeth
TC27-195 15 cm
 curved serrated



Iris Tissue Forceps
TC27-196 10 cm
 1x2 Teeth. Straight



Treves Tissue Forceps
1x2 Teeth
TC27-199 12.5 cm
TC27-200 15.0 cm
TC27-201 18.0 cm



Waugh Tissue Forceps
1x2 Teeth
TC27-202 15 cm
TC27-203 18 cm



Charnley Tissue Forceps
1x2 Teeth
TC27-204 14.5 cm
TC27-205 18.0 cm
 block on jaws.



Dissecting Dressing
Forceps
TC27-206 12.5 cm
TC27-207 15.0 cm
TC27-208 18.0 cm



Dissecting Dressing
Forceps
TC27-209 12.5 cm
TC27-210 15.0 cm
TC27-211 18.0 cm
TC27-212 20.0 cm
 block end, serrated
 jaws.

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط




Ferris-Smith Tissue Forceps 1x2 Teeth
 TC27-213 12 cm
 TC27-214 14 cm
 TC27-215 16 cm
 TC27-216 18 cm
 TC27-217 20 cm




Ferris-Smith Tissue Forceps 2x3 Teeth
 TC27-218 12 cm
 TC27-219 14 cm
 TC27-220 16 cm
 TC27-221 18 cm
 TC27-222 20 cm



Bonney Tissue Forceps
 TC27-223 17 cm



TC27-224
 18 cm
 1x2 Teeth



TC27-225
 18 cm
 2x3 Teeth



Durante Tissue Grasping Forceps
 TC27-226 14.5cm
 TC27-227 16 cm
 TC27-228 20 cm




Russian-Modell Tissue Forceps
 TC27-229 15 cm
 TC27-230 20 cm
 TC27-231 25 cm



Singley-Tuttle Forceps
 TC27-232 23 cm



Collin-Duval Forceps




Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقط



Lucae

TC 27-236
14 cm
Dressing, Serrated
TC 27-237
16 cm
Dressing, Serrated
TC 27-238
18 cm
Dressing, Serrated
TC 27-239
20 cm
Dressing, Serrated



Lucae

1x2 teeth
TC 27-240
14 cm
TC 27-241
16 cm
TC 27-242
18 cm
TC 27-243
20 cm



Adson Bayonet
Dressing Forceps
TC27-244 21 cm
serrated



Yasargil
TC27-245
20 cm



Hunt
TC27-246
20 cm



Hunt
TC27-247
20 cm

Пинцеты
Précèlles - Brucelles
Zahnpinzetten

Пинцети
Pinzas para algodón
ملاقط



Jansen
TC27-248
18.5 cm




Jansen
TC 27-249
16 cm



Adson
TC27-250
23 cm



Troeltsch
TC27-251
12 cm



Troeltsch
TC 27-252
11.5 cm
TC 27-253
14 cm
TC 27-254
17 cm



Landolt
TC27-255
21.5 cm

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



2.7 mm

Debakey Atraumatic Dressing Forceps
 TC27-256 15cm
 TC27-257 20cm
 TC27-258 24cm
 TC27-259 30cm



3.5 mm

Debakey Atraumatic Dressing Forceps
 TC27-260 15cm
 TC27-261 20cm
 TC27-262 24cm
 TC27-263 30cm



2 mm

Cooley Atraumatic Dressing Forceps
 TC27-264 15cm
 TC27-265 20cm
 TC27-266 24cm
 TC27-267 30cm

1.2 mm

TC 27-268
 18 cm
 27-269
 21 cm

2 mm

TC27-270
 18 cm
 TC27-271
 21 cm

Debakey Dressing Forceps



Debakey Forceps



TC27-272 16 cm 1.5 mm
 TC27-273 20 cm 1.5 mm
 TC27-274 24 cm 1.5 mm



TC27-279 16 cm 2.7mm
 TC27-280 20 cm 2.7mm
 TC27-281 24 cm 2.7mm
 TC27-282 30 cm 2.7mm




TC27-275 16 cm 2 mm
 TC27-276 20 cm 2 mm
 TC27-277 24 cm 2 mm
 TC27-278 30 cm 2 mm



TC27-283 16 cm 3.5mm
 TC27-284 20 cm 3.5mm
 TC27-285 24 cm 3.5mm
 TC27-286 30 cm 3.5mm




De Bakey Forceps



TC27-287
 16 cm
 Angular



TC27-289
 24 cm
 Angular



TC27-288
 20 cm
 Angular



TC27-290
 30 cm
 Angular

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقف



Feilchenfeld Forceps

TC27-291
7.5 cm
TC27-292
9 cm
TC27-293
11 cm
Straight



Feilchenfeld

TC27-294
7.5 cm

TC27-295
9.5 cm

TC 27-296
11.5 cm



Hunter Forceps

TC27-297
10.5 cm
Straight



Hunter Forceps

TC27-298
10.5 cm
Curved



Carmalt Forceps

TC27-299
10.5 cm
Straight




Carmalt Forceps

TC27-300
10.5 cm
Curved

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten

 Пінцети
 Pinzas para algodón
 ملاقط



Micro Surgery Instrument



Пинцеты для удаления осколков

Pinces à souder - Pincas et pincas à échardes

Löt-Zange - Pinzette & SplinterZange

Пінцети для видалення осколків

Alicates Soldadura - Pinzas y pinzas de Splinter

لحام كمامية - ملاقط منشقة الملائق



TC 27-311
Soldering Pliers
16cm



TC 27-312
Soldering Pliers
16cm



TC 27-313
Posigrip Articulating
Paper Forceps
15cm



TC 27-314
Standard Articulating
Paper Forceps, Economy Grade
15cm



TC 27-315
Miller
15cm
Curved



Beer Cilia Forceps
TC 27-316
9 cm



Kuehne Cilia Forceps
TC 27-317
10.5 cm



Kuehne Cilia Forceps
TC 27-318
10.5 cm
angled sideways

Пинцеты для удаления осколков


Pinces à souder - Pincas et pincas à échardes

Löt-Zange - Pinzette & SplinterZange

Пінцети для видалення осколків

Alicates Soldadura - Pinzas y pinzas de Splinter

لحام كمامية - ملاقط منشقة الملائق



Splinter Forceps

TC27-319 8 cm
TC27-320 9 cm
TC27-321 10.5 cm
TC27-322 11.5 cm
TC27-323 12.5 cm
TC27-324 14.5 cm

Straight serrated with guide pin

Splinter Forceps

TC27-325 8 cm
TC27-326 9 cm
TC27-327 10.5 cm
TC27-328 11.5 cm
TC27-329 12.5 cm
TC27-330 14.5 cm

Curved serrated with guide pin

Splinter Forceps

TC27-331 8 cm
TC27-332 9 cm
TC27-333 10.5 cm
TC27-334 11.5 cm
TC27-335 12.5 cm
TC27-336 14.5 cm

Straight plain with guide pin

Splinter Forceps

TC27-337 8 cm
TC27-338 9 cm
TC27-339 10.5 cm
TC27-340 11.5 cm
TC27-341 12.5 cm
TC27-342 14.5 cm


Curved plain with guide pin

Stieglitz Forceps with Magnifying Mirrors

TC27-343
12.5 cm




Stieglitz Splinter Forceps
TC27-344 14 cm Str.
TC27-345 14 cm 45° Angled



Stieglitz Splinter Forceps
TC27-346
14 cm 90° Angled



Peet Splinter Forceps
TC27-347
11.5 cm



Ralk Splinter Forceps
TC27-348 15 cm Str.
TC27-349 15 cm Cvd.

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقط



Micro Adson



Adson



Adson



Micro Adson
Light Weight



Adson
Light Weight



Adson
Light Weight



Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقظ



Micro Adson Tissue Forceps
1x2 teeth



TC27-360
12 cm



TC27-361
15 cm



Adson Tissue Forceps
1x2 teeth



TC27-362
12 cm



TC27-363
15 cm



Adson Brown Forceps
3x4 teeth



TC27-364
12 cm



Micro Adson Light Weight
1x2 teeth



TC27-365
12 cm



TC27-366
15 cm



Adson Light Weight
1x2 teeth



TC27-367
12 cm



TC27-368
15 cm




Adson-Brown Light Weight
3x4 teeth



TC27-369
12 cm

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пінцети
Pinzas para algodón
ملاقط



T.C Standard Dressing Forceps



T.C Cushing Forceps
TC27-372 18 cm



T.C Cushing Dressing Forceps
TC27-373
18 cm Cross Serrated



T.C Cushing Dressing Forceps
TC 27-374 18 cm



T.C Dissecting Tissue Forceps
TC27-376 18 cm



T.C Resano Dressing Forceps
TC27-377 25 cm



T.C Cushing Dressing Forceps
TC27-378 19 cm

 Пинцеты
 Précelles - Brucelles
 Zahnpinzetten


 Пінцети
 Pinzas para algodón
 ملاقط



Wangenstein Dressing Forceps T.C.




TC27-379
15.5 cm



TC27-380
23 cm



TC27-381
25.5 cm




T.C DeBakey Micro Needle Tissue Forceps
TC 27-382 18 cm



T.C Bonney Dissecting & Tissue Forceps



TC27-383
1x2 Teeth
18 cm
Straight
X-Serrations



TC27-384
1x2 Teeth
18 cm
Straight, Serrated



TC27-385
2x3 Teeth
18 cm
Straight, Serrated



T.C POTTS-SMITH Dressing Forceps
TC 27-386 15 cm
TC 27-387 18 cm
TC 27-388 20 cm
TC 27-389 23 cm
TC 27-390 25 cm



T.C POTTS-SMITH Dressing Forceps
TC 27-391 16 cm
TC 27-392 18 cm
TC 27-393 21 cm
TC 27-394 24 cm
TC 27-395 26 cm



T.C Evans Pelvic Needle Dressing Forceps
TC 27-396 25.5 cm



T.C Evans Dressing Forceps
TC 27-397 25 cm

Пинцеты
Précettes - Brucelles
Zahnpinzetten

Пинцети
Pinzas para algodón
ملاقط



T.C Mc Indoe Dressing Forceps
TC 27-398 15 cm



T.C De Bakey Dressing Forceps
TC 27-399 15 cm
TC 27-400 19 cm
TC 27-401 24 cm



T.C Gillies Tissue Forceps
TC 27-402 15 cm



Semken
with dissector end

TC 27-403
15 cm



Semken
1x2 teeth

TC 27-404
16 cm

TC 27-405
18 cm



T.C Oehler Tissue Forceps

TC 27-406
14 cm

TC 27-409
20 cm

TC 27-407
16 cm

TC 27-410
25 cm

TC 27-408
18 cm


TC 27-411
30 cm

Комплект для обследования
Ensemble d'examen
Prüfungssatz

Екзаменаційний набір
Conjunto de examen
مجموعة الفحص

TC27-412

First Aid Examination kit Set of 3 Pieces



TC27-413

First Aid Examination Set of 4 Pieces



TC27-414

First Aid Examination Set of 3 Pieces




TC27-415

First Aid Examination Set of 3 Pieces




TC27-416

First Aid Examination Set of 3 Pieces



TC27-417

First Aid Examination Set of 3 Pieces




Комплект для обследования
Ensemble d'examen
Prüfungssatz

Екзаменаційний набір
Conjunto de examen
مجموعه الفحص

TC27-418

First Aid Examination kit Set of 3 Pieces



TC27-419

First Aid Examination Set of 4 Pieces




TC27-420

First Aid Examination Set of 3 Pieces



TC 27-421

First Aid Examination Set of 3 Pieces



TC 27-422

First Aid Examination Set of 3 Pieces



TC27-423

First Aid Examination Set of 3 Pieces




Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين

Single Ended Probe Solid Handle Ø 6 mm, Length 150mm




Fig. 2
Solid Handle



Single Ended Probe Hollow Handle Ø 8 mm Length 150mm



Fig. 2
Anatomic Handle




Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين

EXPLORERS INSTRUMENTS AVAILABLE IN
Hollow Handle Ø 8 mm, Ø 10mm Solid Handle Ø 6mm
Length 170MM

Double-Ended Explorers


EXPLORERS



Зонд
Sonde
Sonde
Sonde

Зонд
Sonda
Sonda
مستقيار


Available in Hollow Handle
Ø 8 mm, Ø 10mm Solid Handle Ø 6mm



Зонды
 Sondes
 Zahnsonden

Зонди
 Sondas
 المستكشفين


Explorers Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




Зонды
 Sondes
 Zahnsonden

Зонди
 Sondas
 المستكشفين

Pointed Probes Available in Hollow Handle Ø 8 mm. Ø 10 mm Solid Handle Ø 6 mm



Pointed Probes Available in Hollow Handle Ø 8 mm. Ø 10 mm Solid Handle Ø 6 mm




Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

BM27-490

Periodontal Set of 10 Pieces Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



TC27-491
Mouth Mirror

TC 27-492
EXP,23/17-A

TC27-493
Exavators

TC27-494
Filling
Instruments

TC27-495
EXP,23/CP-18

TC 27-496
EXP,23/CP-QW

TC 27-497
EXP,23/CP-11

TC 27-498
EXP,23/CP-10

TC27-499
EXP,23/CP-12

TC27-500
EXP,23/CP-11.5-B

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

TC27-501
Columbia Set

TC27-505
McCall Set

Columbia 2L/2R

Columbia 4L/4R

Columbia 13/14

McCalls 14S/18S

McCalls 17/18

McCalls 13S/14S

McCalls 13/14

TC 27-502

TC 27-503

TC 27-504

TC 27-506

TC 27-507

TC 27-508

TC27-509


Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين

TC27-510

William Probe Set Double Ended With Hollow Handle Ø 10 Mm, Length 170mm

Plasma Coated Handle



TC27-511
EXP, 23/CO,QW

TC27-512
EXP, 23/CP,10

TC27-513
EXP, 23/CP,18

TC27-514
EXP, 23/CP,12

TC 27-515
EXP, 23/CP,2

TC27-516
EXP, 23/UNC,15


TC27-517
EXP, 23/CP,11

TC27-518
EXP, 23/CP, 12 .S

Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين

Explorers Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
PLUGGERS




Зонды пародонтологические зонди пародонтологічні
 Sondes parodontales Sondas parodontales
 sonda parodontali المعيار اللثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SPREADERS

HEAT CARRIERS




Зонды пародонтологические
Sondes parodontales
sonde parodontali

зонди пародонтологічні
Sondas parodontales
المسبار اللثة

Endodontics Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm


ANTERIOR SINGLE ENDED



Зонды пародонтологические зонди пародонтологічні
 Sondes parodontales Sondas parodontales
 sonde parodontali المسبار اللثة


Endodontics Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
 POSTERIOR SINGLE ENDED

RETRO FILLING PLUGGERS




■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



■ Инструменты для реставрации
 ■ Instruments pour plombage
 ■ Füllungsinstrumente (mit Kontrawinkel)
 ■ Инструменти для реставрації
 ■ Instrumentos para obturaciones
 ■ ملء الصكوك

Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plomage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Amalgam and Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Інструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Amalgam Instruments and Burnishers Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Інструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Amalgam Instruments and Burnishers Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Plastic Filling Instruments Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Composite Gold Coated Instruments



TC27-640
Markley 2



TC27-641
Westcott N.21



TC27-642
NYU/B2



TC27-643
Filling instruments



TC27-644
Ball Filling Instruments



TC27-645
Filling Instruments



TC27-646
Filling Instruments



TC27-647
Filling Instruments



TC27-648
Plastic Filling Instruments



TC27-649
Heidemann



TC27-650
Plastic Filling Instruments

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Composite Gold Tip Instruments




Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Heidemann Spatula

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Excavators

Excavators Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Composite Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты пародонтологические

Instruments pour la parodontologie

Instrumente für die Parodontologie

Інструменти пародонтологічні

Instrumentos para periodoncia

الصكوك اللثة

Crane Kaplan

Available in Hollow Handle Ø 8 mm



TC 27-690



TC 27-691
Jaquette



TC 27-692
Columbia 13/14



TC 27-693
Kirkland 15/16



TC 27-694
Goldman Fox




TC 27-695
Gracey 11/12



TC 27-696
Gracey 13/14



TC 27-697
Columbia 2R/2L




TC 27-698
Urban Knife 1/2
Lucas Currettes

Available in Hollow Handle Ø 8 mm



TC 27-699
Ø 2 mm



TC 27-700



TC 27-701
Ø 2.5 mm



TC 27-702
Ø 3 mm

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Sinus Lift Instruments Blue Plasma Coating

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



TC27-703
Fig. 1



TC27-704
Fig. 2



TC27-705
Fig. 3




TC27-706
Fig. 4



TC27-707
Fig. 5



TC27-708
Fig. 6




TC27-709
Fig. 7



TC27-710
Fig. 8

Gingival Retractors

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



TC27-711



TC27-712




TC27-713

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Blue Plasma Coating & Operative instruments Available in Hollow Handle Ø 8 mm



Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)


Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

CARVERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plomage

Füllungsinstrumente (mit Kontrawinkel)


Інструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

CARVERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації


Instrumentos para obturaciones

ملء الصكوك

CARVERS INSTRUMENTS


Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

CLEIDS



Nystrom

CHISELS



Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Інструменти для реставрації


Instrumentos para obturaciones

ملء الصكوك


ENAMEL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

MARGINTRIMMERS



HATCHES



Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm
CHISELS



Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Інструменти для реставрації


Instrumentos para obturaciones

ملء الصكوك

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

Heidemann Spatulas



SPATULAS



Экскаваторы
Excavateurs
Excavatoren


Εκσκαβатори
Excavadores
حفارات

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm




Экскаваторы
Excavateurs
Excavatori
Excavadores
حفارات

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Oval Shape



Экскаваторы
Excavateurs
Excavatoren

Εκσκαβатори
Excavadores
حفارات

Spoon Excavators Round Shape

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



125/126



TC 27-814



127/128



TC 27-815



129/130



TC 27-816



133/134



TC 27-817



153/154



TC 27-818



3 NYU E1



TC 27-819



4 NYU E2



TC 27-820



5 NYU E3




TC 27-821

Зонды
Sondes
Zahnsonden

Зонди
Sondas
المستكشفين

PROBES



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة



TC 27-830
Excavator



TC 27-831
Excavator



TC 27-832
Excavator




TC 27-833
Filling Instruments




TC 27-834
Filling Instruments



TC 27-835
Filling Instruments



TC 27-836
Filling Instruments



TC 27-837
Filling Instruments



TC 27-838
Cement Spatula

Инструменты пародонтологические

Instruments pour la parodontologie




Instrumente für die Parodontologie

Інструменти пародонтологічні

Instrumentos para periodoncia

الصكوك اللثة

MODELLING INSTRUMENTS



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Periodontal Instruments Set



TC27-847
Set of 10 Pieces.

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة


Periodontal Files

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

BM27-848

Sugarman File Set

Excellent to remove and to crush strong calculus deposits



TC 27-849
9/10
Schluger
Curved File

TC 27-850
11/12
Buck
Straight File

TC 27-851
1S/2S
Sugarman File

TC 27-852
3S/4S
Sugarman File

TC 27-853
3-7
Hirschfeld

TC 27-854
5-11
Hirschfeld

TC 27-855
9-10
Hirschfeld

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة


Periodontal Files

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

BM27-856

Sugarman File Set

Excellent to remove and to crush strong calculus deposits



TC27-857
9/10
Schluger
Curved File

TC27-858
11/12
Buck
Straight File

TC27-859
1S/2S
Sugarman File

TC27-860
3S/4S
Sugarman File

TC27-861
3-7
Hirschfeld

TC27-862
5-11
Hirschfeld


TC 27-863
9-10
Hirschfeld

Кюрети
Curetas
Küretten

Кюрети
Curette
كشش

SICKLE SCALERS

Designed for the removal of supragingival calculus



TC27-864
H6/H7
Hygienist

TC27-865
H5/33
Hygienist-
Jacquette

TC27-866
H5/U15
Hygienist-
Towner

TC27-867
6/7
Goldman

TC27-868
U15
Towner

TC27-869
U15/30
Towner-
Jacquette


TC27-870
U15/33
Towner-
Jacquette

Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین

SCALERS

JACQUETTE SCALERS



TC27-871
Nevi 1
Anterior

TC27-872
Nevi 1/H5
Anterior

TC 27-873
Nevi 2
Posterior

TC27-874
Nevi 3
Posterior

TC27-875
Nevi 4
Posterior

TC 27-876
30/33
Jacquette

TC 27-877
31/32
Jacquette


TC27-878
34/35
Jacquette

TC27-879
H3/H4
Jacquette

TC 27-880
14/15
Posterior

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتیس



Кюрети
Curetas
Küretten

Кюрети
Curette
کرتیس

SICKLE SCALERS

POSTERIOR

ANTERIOR



TC 27-890
135

TC 27-891
REM 3/4
Remington

TC 27-892
129

TC 27-893
T2/T3
Taylor

TC 27-894
137

TC 27-895
0/00
Morse

Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین



TC 27-896
137
Mini Five



TC 27-897
OD1/33
Odont/
Jacquette



TC 27-898
128
Nebraska



TC 27-899
USC-128



TC 27-900
WH-2
Whiteside




TC 27-901
YG-15
Younger-
Good

Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین

BALL SCALERS

CRANE-KAPLAN



TC 27-902
CI-1

TC27-903
CI-2/3

TC 27-904
152
Cumine

TC27-905
B-1

TC27-906
B-4

TC27-907
B-2/3

TC27-908
B-5/6

TC27-909
CK-6

TC27-910
CK-6D


TC27-911
333/334
Mini Kaplan

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتیس

SCALER / CURETTE SETS

Mc CALL SET



TC27-912
11/12
Scaler

TC27-913
11A/12A
Scaler

TC27-914
11/12S
Modified
McCall

TC 27-915
13S/14S
Pointed
Curette

TC27-916
MTC 13/14
Pointed
Curette

TC 27-917
17/18
Curette

TC 27-918
17S/18PT
Curette

TC 27-919
19/20
Curette


TC27-920
Goldman-Fox

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتیس

SCALER / CURETTE SETS

GOLDMAN-FOX SET



TC27-921
GF-1
Universal Scaler

TC27-922
GF-2
Curette

TC27-923
GF-3
Curette

TC27-924
GF-4
Curette

TC27-925
GF-5 Hoe
Buccal/Lingual

TC27-926
GF-6
Hoe

TC27-927
GF-21

TC27-928
GF-21S/33
Combination

Кюрети
Curetas
Küretten


Кюрети
Curette
کرسین

SCALER / CURETTE SETS

GOLDMAN-FOX SET

CATTONI SET

COLUMBIA



Кюрети
Curetas
Küretten

Кюрети
Curette
كشش

UNIVERSAL CURETTES

BARNHART



TC 27-939
IU13/14 Indiana
University
Anterior



TC 27-940
IU17/18 Indiana
University
Posterior



TC 27-941
1/2 Bunting
Curette-Chisel
combination



TC 27-942
5/6 Bunting



TC 27-943
YG7/8
Younger-Good



TC 27-944
1/2 Barnhart
(shown)

TC 27-945
1/2 Barnhart
Rigid



TC 27-946
5/6 Barnhart
(shown)

TC 27-947
5/6 Barnhart
Rigid

Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین



TC27-948 1/2
(shown)

TC27-949 1/2
Rigid



TC27-950 3/4
(shown)

TC27-951 3/4
Rigid



TC27-952 5/6
(shown)

TC27-953 5/6
Rigid

TC27-954 5/6
Extra Rigid



TC27-955 7/8
(shown)

TC27-956 7/8
Rigid

TC 27-957 7/8
Extra Rigid



GRACEY 1/2



GRACEY 7/8



GRACEY 11/12



GRACEY 13/14



GRACEY 15/16



GRACEY 17/18



TC27-958 9/10
(shown)

TC27-959 9/10
Rigid



TC27-960 11/12
(shown)

TC27-961 11/12
Rigid

TC27-962 11/12
Extra Rigid



TC27-963 13/14
(shown)

TC27-964 13/14
Rigid

TC27-965 13/14
Extra Rigid



TC 27-966 15/16
(shown)

TC27-967 15/16
Rigid



TC27-968 17/18
(shown)

TC27-969 17/18
Rigid


Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین

GRACEY CURETTES

MICRO MINI

GRACEY PROPHY CURETTES



TC27-970
Micro Mini
Five 1/2

TC 27-971
Micro Mini
Five 7/8

TC27-972
Micro Mini
Five 11/12

TC27-973
Micro Mini
Five 13/14

TC27-974
P1/2

TC27-975
P3/4


TC27-976
P5/6

TC27-977
P7/8

Кюрети
Curetas
Küretten

Кюрети
Curette
كشش

Gracey Curettes Solid Handle Ø 6 mm, Length 170mm



TC27-978
1/2

TC 27-979
3/4

TC 27-980
5/6

TC27-981
7/8

TC27-982
9/10

TC 27-983
11/12

TC 27-984
13/14


TC 27-985
15/16

TC27-986
17/18

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتن

Gracey Curettes Hollow Handle Ø 10 mm, Length 170mm



Кюрети
 Curetas
 Küretten

Кюрети
 Curette
 كرشين



Hollow Handle Ø 10 mm



Hollow Handle Ø 10 mm

TC27-1003
Gracey Curette Set of 7 Pcs
Gracey Curette 1-2
Gracey Curette 3-4
Gracey Curette 5-6
Gracey Curette 7-8
Gracey Curette 9-10
Gracey Curette 11-12
Gracey Curette 13-14

TC27-1011
Gracey Curette
(Black & Gold Coating)
Set of 7 Pcs
Gracey Curette 1-2
Gracey Curette 3-4
Gracey Curette 5-6
Gracey Curette 7-8
Gracey Curette 9-10
Gracey Curette 11-12
Gracey Curette 13-14

- Кюрети
- Curetas
- Küretten




TC27-1027
Gracey Curette Set of 7 Pcs
with Silicon Handle
Gracey Curette 1-2
Gracey Curette 3-4
Gracey Curette 5-6
Gracey Curette 7-8
Gracey Curette 9-10
Gracey Curette 11-12
Gracey Curette 13-14

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتس

HOLLOW HANDLE PLASMA BLUE COATED Ø 10 MM, LENGTH 170MM



TC27-1035

Gracey Curettes Set of 7 Pcs

Gracey Curettes 1/2
Gracey Curettes 3/4
Gracey Curettes 5/6
Gracey Curettes 7/8
Gracey Curettes 9/10
Gracey Curettes 11/12
Gracey Curettes 13/14

Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین

HOLLOW HANDLE Ø 10 MM, LENGTH 170MM WITH SILICON RING



TC27-1043

Gracey Curettes Set of 7 Pcs

Gracey Curettes 1/2

Gracey Curettes 3/4

Gracey Curettes 5/6

Gracey Curettes 7/8

Gracey Curettes 9/10

Gracey Curettes 11/12

Gracey Curettes 13/14

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتش



TC27-1044
1/2



TC27-1045
3/4



TC27-1046
5/6



TC27-1047
7/8



TC27-1048
9/10



TC27-1049
11/12



TC27-1050
13/14

TC27-1051

**Gracey Curettes Set of 7 Pcs
with Silicon Grip**

Gracey Curettes 1/2

Gracey Curettes 3/4

Gracey Curettes 5/6

Gracey Curettes 7/8

Gracey Curettes 9/10


Gracey Curettes 11/12

Gracey Curettes 13/14

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتین

Silicone Handles



TC 27-1052
Explorer 1

TC 27-1053
Explorer 2

TC 27-1054
Explorer 3

TC 27-1055
Root Canal Spreader 1

TC 27-1056
Root Canal Spreader 2

TC 27-1057
Root Canal Spreader 3

TC 27-1058
Root Canal Spreader 4



TC 27-1059
Amalgam
Instruments 1



TC 27-1060
Amalgam
Instruments 2



TC 27-1061
Amalgam
Instruments 3



TC 27-1062
Amalgam
Instruments 4



TC 27-1063
Amalgam
Instruments 5



TC 27-1064
Amalgam
Instruments 6

Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین



TC27-1065
Filling Instruments 1



TC27-1066
Filling Instruments 2



TC27-1067
Filling Instruments 3



TC27-1068
Filling Instruments 4



TC27-1069
Filling Instruments 5



TC27-1070
Lucas Curettes 1



TC27-1071
Lucas Curettes 2



TC27-1072
Lucas Curettes 3



TC27-1073
Lucas Curettes 4




TC27-1074
Lucas Curettes 5

Кюрети
Curetas
Küretten

Кюрети
Curette
کرتس

GRACEY MINI FIVE
 HOLLOW HANDLE Ø 10 MM, LENGTH 170MM
 CERAMIC COATED HANDLE

GRACEY AFTER FIVE
 HOLLOW HANDLE Ø 10 MM, LENGTH 170MM
 CERAMIC COATED HANDLE



TC27-1075
 Gracey
1/2

TC27-1076
 Gracey
3/4

TC27-1077
 Gracey
5/6

TC27-1078
 Gracey
7/8

TC27-1079
 Gracey
11/12

TC27-1080
 Gracey
13/14

TC27-1081
 Gracey
1/2 AF

TC27-1082
 Gracey
3/4 AF

TC27-1083
 Gracey
5/6 AF

TC27-1084
 Gracey
7/8 AF

TC27-1085
 Gracey
11/12 AF


TC27-1086
 Gracey
13/14 AF

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Sickle Scalers Instruments

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

JACQUETTE SCALERS

CRANE-KAPLAN



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SCALER CURETTE



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

SICKLE SCALERS INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

SICKLE SCALERS



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie


Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

PERIODONTAL INSTRUMENTS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

JACQUETTE SCALERS

CRANE-KAPLAN



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

PERIODONTAL INSTRUMENTS
Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm


SCALER CURETTE



Стоматологический набор скейлеров
Instruments a detartre
Zahnreiniger

Стоматологічний набір скейлерів
Limpiadores de sarro
المستظفون


SCALERS Solid Handle Ø 6 mm



Стоматологический набор скейлеров
Instruments a detartre
Zahnreiniger

Стоматологічний набір скейлерів
Limpiadores de sarro
المستظفون

Tartar Remover Instruments Solid Handle Ø 6 mm




Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

PERIODONTAL KNIVES



TC27-1142
1
Reddy Thin

TC27-1143
2
Reddy Thin

TC27-1144
3/4
Buck

TC27-1145
5/6
Buck

TC27-1146
1/2
Solt

TC27-1147
1/2
Sanders

TC27-1148
USC Towner
19/20

TC27-1149
CK-3
Crane Kaplan


TC27-1150
Dr. Allen
End-Cutting
Intrasulcular Knife

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

PERIODONTAL KNIVES



TC27-1151
15K/16K
Kirkland

TC27-1152
1/2
Orban

TC27-1153
1/2 Orban
Knife

TC27-1154
7
Goldman-Fox

TC27-1155
8
Goldman-Fox

TC27-1156
9
Goldman-Fox

TC27-1157
11
Goldman-Fox

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie


Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

PERIODONTAL CHISELS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

OCHSENBEIN

FEDI



TC27-1158
1
Ochsenbein

TC27-1159
2
Ochsenbein

TC27-1160
3
Ochsenbein

TC27-1161
4
Ochsenbein

TC27-1162
1
Fedi

TC27-1163
2
Fedi

TC27-1164
3
Fedi


TC27-1165
4
Fedi

инструменты для костной ткани
Limes a os
Knochenfeilen

інструменти для кісткової тканини
Limas para huesos
ملفات العظام



TC27-1166
Fig. 1
PoloKoff



TC27-1167
Fig. 33



TC27-1168
Fig. 45
Miller



TC27-1169
Fig. 64
Miller



TC27-1170
Fig. 10



TC27-1171
Fig. 1
Wahl



TC27-1172
Fig. 2
Wahl



TC27-1173
Fig. 12
Howard



TC27-1174
Fig. 12A




TC27-1175
Fig. 12CA

Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوك اللثة

PERIODONTAL CHISELS

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



TC27-1176
1/2
Kramer-Nevins

TC27-1177
36/37
Rhodes
Back-Action

TC27-1178
44/45
Remington
Hoe

TC27-1179
TG

TC27-1180
TG-0
Modified TG with
Ochsendein cutouts

TC27-1181
3/4
Solt

TC27-1182
13K/TG

TC27-1183
13K/13KL

Кюрети
Curetas
Küretten


Кюрети
Curette
کرسین

BONE CURETTES

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

OVAL

ROUND



TC27-1184
Hemingway

TC27-1185
Hemingway

TC27-1186
Hemingway

TC27-1187
Martin

TC27-1188
Martin

TC27-1189
Martin


Кюрети
Curetas
Küretten

Кюрети
Curette
کرسین

LUCAS SPOON SHAPE

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

LUCAS



TC27-1190
84

TC27-1191
85

TC27-1192
86

TC27-1193
87

TC27-1194
88

TC27-1195
85/86
Lucas
Surgical Curette

TC27-1196
86/87
Lucas
Surgical Curette

TC27-1197
Lucas

TC27-1198
Lucas

Инструменты для реставрации

Instruments pour plombage

Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

P.K. THOMAS WAXING INSTRUMENTS

P.K. THOMAS WAX CARVER SET OF 5 PCS


Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm

P.K. THOMAS WAX CARVER SET OF 5 PCS

Aluminum Anodizing Colours Handle Ø 8 mm



TC27-1204
P.k. Thomas Wax Carvers Set of 5 Pcs



TC27-1210
P.k. Thomas Wax Carvers Set of 5 Pcs
with Sterilized Cassette

Инструменты для реставрации

Instruments pour plombage


Füllungsinstrumente (mit Kontrawinkel)

Инструменти для реставрації

Instrumentos para obturaciones

ملء الصكوك

COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC



TC27-1211

TC27-1212

TC27-1213

TC27-1214

TC27-1215

TC27-1216



TC27-1217
Composite Instruments Blue
Titanium with Cassettes

Шпатели для цемента
Spatules à ciment
Zementspatel

Шпателі для цементу
Espátulas para cemento
ملاعق الاسمنت

Available in Hollow Handle Ø 8 mm, Ø 10 mm Solid Handle Ø 6 mm



TC27-1218
Fig. 1



TC27-1219
Fig. 2



TC27-1220
Fig. 3



TC27-1221
Fig. 00



TC27-1222
Fig. 0



TC27-1223
Fig. 1



TC27-1224
Fig. 2



TC27-1225
Fig. 3



TC27-1226
Fig. 4



TC27-1227
Weston



TC27-1228




TC27-1229

Шпатели для цемента
Spatules à ciment
Zementspatel

Шпателі для цементу
Espátulas para cemento
ملاعق الاسمنت

Wax and Modelling Carvers



TC27-1230
Le Cron

TC27-1231
Zahle

TC27-1232
Zahle

TC27-1233
Zahle

TC27-1234
Gritman

TC27-1235
Beale

TC27-1236
Roach

TC27-1237
Vehe

Инструменты для реставрации
Instruments a modeler
Modellierinstrumente

Инструменти для реставрації
Instrumentos para modelar
ابوات تقطیع



BM27-1238

Graded Handle Hylin Carver



TC27-1239
Hylin



TC27-1240
Hylin-Le Cron



TC27-1241
Hylin-Ceramica



TC27-1242
Hylin



TC27-1243



TC27-1244
Le Cron



TC27-1245
Le Cron



TC27-1246
Le Cron

Ножи для воска
Couteaux a cire
Wachsmesser

Ножі для воску
Cuchillos para cera
السكاكين الشمع

WAX KNIVES



TC27-1247
Fig. 1
Fahnenstock
125mm

TC27-1248
Fig. 2
Fahnenstock
170mm

TC27-1249
Fig. 1
Lessmann
125mm

TC27-1250
Fig. 2
Lessmann
170mm

TC27-1251
Fig. 1
with small spoon
125mm

TC27-1252
Fig. 2
with large spoon
170mm

Ножи для воска
Couteaux a cire
Wachsmesser

Ножі для воску
Cuchillos para cera
السكاكين الشمع

PLASTER KNIVES



TC27-1253
160mm



TC27-1254
160mm



TC27-1255
Buffalo



TC27-1256
Buffalo



TC27-1257
Plaster Knife



TC27-1258
170mm

Ножи для воска
Couteaux a cire
Wachsmesser

Ножі для воску
Cuchillos para cera
السكاكين الشمع

WAX KNIVES



TC27-1259
Sachs



TC27-1260
130mm
Ermert



TC27-1261
Hook Knives



TC27-1262
125mm
Antkoviak



BM27-1263



Special Inox
TC27-1264 170mm
TC27-1265 200mm

Шпатель для замешивания
Spatules a alginate et a platre
Alginat- und Gipsspatel

Шпатель для змішування
Espatulas para alginato y para yeso
ملاعق الجبس الجينات

WAX SPATULAS



TC27-1266
Fig. 1
18.5 cm



TC27-1267
Fig. 2
18 cm



TC27-1268
Wooden Handle
19.5cm curved



TC27-1269
Wooden Handle
20cm curved



TC27-1270
Wooden Handle
22.5cm curved





TC27-1271
Wooden Handle
19.5cm



TC27-1272
Wooden Handle
21.5cm





TC27-1273
Wooden Handle
20.5cm



Compressor Single Flask
TC27-1278


Compressor Double Flask
TC27-1279

Compressor Double Flask
TC27-1280



Артикуляторы
Artocilar
Artikulations

Артикулятори
Articuler
ارتیکولاتورس



TC27-1286
Articulator



TC27-1287
Articulator




TC27-1288
Articulator
S



TC27-1289
Articulator
M




TC27-1290
Articulator
L



TC27-1291
Articulator

Артикуляторы
Artocilar
Artikulations

Артикулятори
Articuler
ارتیکولیٹوری



TC27-1292
Articulator



TC27-1293
Articulator



TC27-1294
Articulator



TC27-1295
Articulator



TC27-1296
Articulator



TC27-1297
Articulator

■ Артикуляторы
■ Artocilar
■ Artikulations

■ Артикулятори
■ Articulier
■ آرتیکولاتورس


TC27-1298

Technical Datas:

Intercondylar distance 120 mm.
condylar inclination fixed at 30°
Bennet angle fixed at 7°
Incisal table inclination 10°
Distance between mounting plates 110 mm.
Weight 700 gr.
Opening 180°

features:

Lightweight anathomic articulator made of aluminum die-casting. it can be used following the geometrical laws of medium values (Bonwill triangle). Fixing according to centric. Automatic return to the hinge axis with springs. Perfect stability also in rec'versed position. Easy access and visibility to the lingual surfaces. interchangeable self centring aluminum mounting plates.



TC27-1299

Technical Datas:

Intercondylar distance 120 mm.
condylar inclination fixed at 30°
Bennet angle fixed at 7°
Incisal table inclination flat
Distance between mounting plates 110 mm.
Weight 700 gr.
Opening 180°

Special Features:

Adjustable back by registration screw to work without incisal rod by keeping constant trajectories and inclinations during side and back/forward movements. Pure translation movement (key and instructions included). Anti-shock hinge rod in case articulator falls down. Bigger size lock screws for mounting plates.



Increased height and width

- (1) More space when using removable stumps systems.
- (2) More space for installation of mobile prosthesis with high chewing or big base.
- (3) More space when using magnetic plates



TC27-1300

Technical Datas:

Intercondylar distance 120 mm.
condylar inclination fixed at 30°
Bennet angle fixed at 7°
Incisal table inclination 10°
Distance between mounting plates 85 mm.
Weight 400 gr.
Opening (item 5000) 125°
Opening (item 5010) 180°




TC27-1301

Technical Datas:

Intercondylar distance 110 mm.
condylar inclination fixed at 30°
Bennet angle fixed at 7°
Incisal table inclination flat
Distance between mounting plates 85 mm.
Weight 450 gr.
Opening 180°



TC27-1302
Anvil
Dental Laboratory
Instrument




TC27-1303 Adult Palatal /
TC27-1304 Child Palatal

TC27-1305 Buccal

TC27-1306 Adult Palatal /
TC27-1307 Child Palatal

TC27-1308 Adult Palatal /
TC27-1309 Adult Buccal




TC27-1310 Lingual upper /
TC27-1311 Lingual Lower

TC27-1312 Adult Buccal /
TC27-1313 Child Buccal

TC27-1314 Child Palatal /
TC27-1315 Child Buccal

TC27-1316
Palatal-Adult



TC27-1317
Palatal-Adult Large

TC27-1318
Palatal-Child

TC27-1319
Buccal-Small

TC27-1320
Buccal

TC27-1321
Lingual

TC27-1322
Lingual-Large



Wazirbad Road, Ugoki, Sialkot-51310 PAKISTAN

Phone: +92-52-3561347-3553588

Phone: +92-52-3561202

Cell: +92-300-6139615

Cell: +92-322-6139615

Email: info@tcmedicline.com

Email: exporttcmedicline@gmail.com