


MAXILLOFACIAL SURGERY INSTRUMENTS



MANUFACTURER & EXPORTER

MAXIOFACIAL SURGERY INSTRUMENTS

MAXIOFACIAL SUARGERY INSTRUMENTS

A Bard-Parker handle with a No. 15 blade is the usual scalpel used to make incisions on the gums where a flap can be raised to gain access to tooth.

RETRACTORS

Tongue, cheek and flap retractors are needed to provide adequate access to the tooth and are legion in number.

NEEDLE HOLDERS AND SUTURES

Sutures are used to replace and secure the gum flap in its former position, to limit the size of the wound and therefore restrict subsequent contamination with food debris, and to assist in controlling bleeding. The usual suture size for intra-oral work is 4/0 and the material may be absorbable or non-absorbable Scissors

Dissecting scissors and suture-cutting scissors are usually used.

ARTERY FORCEPS

These are designed to hold soft tissue without damaging it. College tweezers are used to lift up sutures prior to removing them.

PERIOSTEAL ELEVATORS

A number are available and are similarly used like a normal elevator but periosteal elevators have two heads – one at each end. One end is sharp-tipped to expand the tooth socket while the other has a blunt end to remove granulation tissues or raise a gum flap.

CHISELS AND MALLET


Chisels are more appropriate in the operating theatre and are particularly useful for disto-angular third molars and upper wisdom teeth extraction.

BONE CURETTES

The Mitchells trimmer is probably the most valuable instrument in the category and is used to remove granulation tissue in the tooth socket.

BONE RONGUER AND BONE FILE

Used to smoothen any sharp bone edges after tooth extraction is done.



MANUFACTURER & EXPORTER

Ручки для скальпелей
Manches de bistouris
Skalpellgriffe

Ручки для скальпелів
Mangos de bisturi
مشرط مقابض

Scalpel Handles



TC 28-001
Fig. 3
12.5cm



TC 28-002
Fig. 3
12.5cm



TC 28-003
Fig. 7
16.5cm



TC 28-004
Fig. 5
14.5cm



TC 28-005
Fig. 5A
14.5cm



TC 28-006
6 Way Adjustable
Angled
Scalpel Handle



TC 28-007
Fig. 3L
21cm



TC 28-008
Fig. 4L
21.8cm



TC 28-009
Fig. 4
13.5cm




TC 28-010
Double Handle
No.4 + No.3
16.3cm


Ножи для скальпеля и скальпеля
Blalpel Blade & Scalpel Knives
Scalpel Blade & Scalpel Messer

Ножи скальпельного леза та скальпеля
Cuchillos de Escalpele y Escalpele
مشروط شفرة و مشروط المسكين


STERIL SCALPEL BLADES



SCALPEL KNIVES



■ Скальпельные ножи ■ Ножі скальпеля
 ■ Couteaux à scalpel ■ Cuchillos Scalpel
 ■ Skalpellmesser ■ سكاكين مشرط



Роторасширители
 Ecarteurs - Ouvre-bouches
 Vestibulumabhalter / Mundspanner

Роторозширювачі
 Separadores - Abre-bocas
 الكاشحات - الكاشحات الفم



TC 28-041 10.5cm
TC 28-042 13.5cm
 Weitlaner
 Fig. 1



TC 28-043 10.5cm
TC 28-044 13.5cm
 Weitlaner
 Fig. 2



TC 28-045
 Stainless Steel
 Double
 Ended Cheek
 Retractors




TC 28-046
 Stainless Steel
 Single
 Ended Cheek
 Retractors




TC 28-047 11cm
TC 28-048 13cm
TC 28-049 15cm
 Weitlander



TC 28-050 11cm
TC 28-051 13cm
TC 28-052 15cm
 Jennings



TC 28-053
 Denhard
 12cm



TC 28-054
 Roser-Konig
 16cm



TC 28-055
 Roser-Konig
 19cm

Роторасширители
 Ecarteurs - Ouvre-bouches
 Vestibulumabhalter / Mundspanner

Роторозширювачі
 Separadores - Abre-bocas
 الكاشحات - الكمامات الفم



TC 28-056
HMG
Heister
13cm



TC 28-057
Mouth Gag Molt
11 cm for (Children)



TC 28-058
Mouth Gag Molt
14cm for (Adults)




TC 28-059
14cm

TC 28-060
15cm

Cheek Retractor



TC 28-061
Wangen Holder
15cm, with or without
lip protection



TC 28-062
Wangen Holder
15cm, with
or without lip protection



TC 28-063
Wangen Holder
Wangen holder with dual
function, 17cm



TC 28-064
Retractor with
Ergo Grip
24cm
ergonomic handle



TC 28-065
Cheek Retractor
Retractor Ivory
Width adjustable
& flexible, 20cm

Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
الكامشات

RETRACTORS



TC 28-066
Columbia
14.5 cm



TC 28-067
Sternberg
14 cm
Double Ended



TC 28-068
Sternberg
16 cm



TC 28-069
Sternberg
14 cm
Single Ended



TC 28-070
Oringer Medium
11x8 cm



TC 28-071
Brunings
Cheek Retractor
19cm



TC 28-072
Minnesota
Small
14cm



TC 28-073
Minnesota
Large
16cm



TC 28-074
Cheek Retractor
20cm
Easy to Bend
by Hand



TC 28-075
Cheek Retractor
and Tongue
Pushers Bishop with
Comfortable Margin,
16cm



TC 28-076
Cheek Retractor
Combination
25cm

Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
الكامشات



TC28-077
Cheek Retractor
Long



TC28-078
Cheek Retractor
Long



TC28-079
Cheek Retractor
Long



TC28-080
Unique Retractor
Length 17.5cm
angled 45°



TC28-081
Unique Retractor
Length 17.5cm
angled 90°



TC28-082
Retractor
16cm



TC28-083
Retractor
19cm



TC28-084
Mini-Langenbeck
Retractor
16cm, Hooks 22 mm



TC28-085
Cheek Retractor
Combination
21cm, 2 Pcs



TC28-086
Senn Miller Retractor
16cm Hooks 22 mm



TC28-087
75x75
Austin Tissue
Retractors

8cm
Available with
round tip on
request

13cm
Available with
round tip on
request



TC28-088
75x75
Austin Tissue
Retractors



Joll (Friedman)
TC28-089 15.5 cm
Thyroid- and Vaginal Retractor


Ранорасширители
 Ecarteurs
 Wundhaken

Ранорозширювачі
 Separadores
 الكاشحات



Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
الكامشات



Heiss Retractor
TC 28-106 10 cm
Sharp 4:4
TC 28-107 10 cm
Blunt 4:4

Finsen Retractor
TC 28-108 7 cm
3 x 4 teeth, sharp
TC 28-109 7 cm
3 x 4 teeth, Blunt

Finsen Retractor
TC 28-110 5 cm
2 x 3 teeth sharp/blunt
TC 28-111 7 cm
2 x 3 teeth sharp/blunt

Buckley Retractor
TC 28-112
6 cm



ALM Retractor
TC 28-113 7 cm
4 x 4 teeth sharp/blunt
TC 28-114 10 cm
4 x 4 teeth sharp/blunt

Jansen Retractor
TC 28-115 10 cm
Sharp 3x3
TC 28-116 10 cm
Blunt 3x3

**Markham-Meyerding
Hemilaminectomy Retractor**
TC 28-117 18cm
Side Blade 2.5cm wide

Gelpi-Loktite Retractor


TC 28-118
9 cm
6 mm Sharp
TC 28-119
11 cm
10 mm Sharp
TC 28-120
14 cm
10 mm Sharp
TC 28-121
16 cm
10 mm Sharp
TC 28-122
19cm
10 mm Sharp

Ранорасширители
Ecarteurs
Wundhaken

Ранорозширювачі
Separadores
المنشآت



Mollison Retractor
TC28-123 15 cm
sharp



Milligan Retractor
TC28-124 13 cm
18.5 mm x 20 mm
Blunt, 3x3



Johnson Retractor
TC28-125 12.5cm



Miskimon Retractor
TC28-126 22 cm
33 mm x 38 mm
Sharp, 7x7



Travers Retractor
TC28-127
21 cm Blunt



Righby Retractor
TC28-128 35 x 20 mm
TC28-129 60 x 20 mm
TC28-130 65 x 15 mm
18 cm



TC28-131 31 cm
3x4 Blunt prongs,
13 mm x 14.5 mm
TC28-132 16.5 cm
17 mm x 19 mm
3x4 Sharp prongs,
TC28-133 16.5 cm
17 mm x 19 mm
3x4 Blunt prongs,
TC28-134 20 cm
21 mm x 210 mm
3x4 Blunt prongs,

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



TC 28-135
Micro Iris
9 cm
Straight



TC 28-136
Micro Iris
9 cm
Curved




TC 28-137
Iris
11.5 cm
Curved



Iris
TC 28-138 9cm
TC 28-139 10.5cm
TC 28-140 11.5cm
TC 28-141 12cm
Straight



Iris
TC 28-142 9cm
TC 28-143 10.5cm
TC 28-144 11.5cm
TC 28-145 12cm
Curved



TC 28-146
Angular. Saw edge
11.5 cm



TC 28-147
Angular
11.5 cm

 Ножицы десневые
 Ciseaux à gencives
 Zahnfleischscheren

 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ



TC 28-148
 Goldman-Fox
 13 cm
 Saw edge,
 straight



TC 28-149
 Goldman-Fox
 13 cm
 Saw edge,
 curved




TC 28-150
 Goldman-Fox
 13 cm
 Saw edge,
 double curved




TC 28-151
 La Grange
 11.5 cm
 Saw edge,
 double curved



TC 28-152
 Neumann
 12 cm
 Saw edge,
 double curved



TC 28-153
 Kilner (Ragnell)
 12 cm
 Straight




TC 28-154
 Kilner (Ragnell)
 12 cm
 Curved



TC 28-155
 Joseph
 14 cm
 Straight

 Ножицы десневые
 Ciseaux à gencives
 Zahnfleischscheren


 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ



Joseph
TC28-156 14 cm
Curved



Kelly
TC28-157 16cm
TC28-158 18cm
Straight



Kelly
TC28-159 16cm
TC28-160 18cm
Curved



Locklin
TC28-161 16 cm
Saw edge,
angular



Locklin
TC28-162 12 cm
TC28-163 16 cm
Saw edge,
angular



TC28-164
Dean
17 cm
Saw edge,
angular



TC28-165
Straight
12 cm



TC28-166 11.5cm
TC28-167 13cm
TC28-168 14cm
Straight




TC28-169 11.5cm
TC28-170 13cm
TC28-171 14cm
Curved




Mayo
TC28-172 14 cm
TC28-173 16 cm
TC28-174 18 cm
Straight
TC28-175 14 cm
TC28-176 16 cm
TC28-177 18 cm
Curved

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



TC28-178
Universal
12 cm
Saw edge



Beebee
TC28-179 10cm
TC28-180 12cm
Saw edge, straight



Beebee
TC28-181 10cm
TC28-182 12cm
Saw edge, curved



TC28-183
15 cm



TC28-184
Dental scissors
large 17.5 cm



TC28-185 14cm
TC28-186 18cm
universal
autoclavable up to
143 °C

 Ножицы десневые
 Ciseaux à gencives
 Zahnfleischscheren

 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ



Standard Operating Scissors

Sharp/Sharp Straight

TC 28-187	11.5 cm
TC 28-188	13.0 cm
TC 28-189	14.0 cm
TC 28-190	15.0 cm
TC 28-191	16.0 cm
TC 28-192	17.5 cm
TC 28-193	18.5 cm
TC 28-194	20.0 cm



Standard Operating Scissors

Sharp/Sharp Curved


TC 28-195	11.5 cm
TC 28-196	13.0 cm
TC 28-197	14.0 cm
TC 28-198	15.0 cm
TC 28-199	16.0 cm
TC 28-200	17.5 cm
TC 28-201	18.5 cm
TC 28-202	20.0 cm



Standard Operating Scissors

Blunt/Blunt Straight

TC 28-203	11.5 cm
TC 28-204	13.0 cm
TC 28-205	14.0 cm
TC 28-206	15.0 cm
TC 28-207	16.0 cm
TC 28-208	17.5 cm
TC 28-209	18.5 cm
TC 28-210	20.0 cm



Standard Operating Scissors

Blunt/Blunt Curved

TC 28-211	11.5 cm
TC 28-212	13.0 cm
TC 28-213	14.0 cm
TC 28-214	15.0 cm
TC 28-215	16.0 cm
TC 28-216	17.5 cm
TC 28-217	18.5 cm
TC 28-218	20.0 cm



Standard Operating Scissors

Sharp/Blunt Straight

TC 28-219	11.5 cm
TC 28-220	13.0 cm
TC 28-221	14.0 cm
TC 28-222	15.0 cm
TC 28-223	16.0 cm
TC 28-224	17.5 cm
TC 28-225	18.5 cm
TC 28-226	20.0 cm




Standard Operating Scissors

Sharp/Blunt Curved

TC 28-227	11.5 cm
TC 28-228	13.0 cm
TC 28-229	14.0 cm
TC 28-230	15.0 cm
TC 28-231	16.0 cm
TC 28-232	17.5 cm
TC 28-233	18.5 cm
TC 28-234	20.0 cm

■ Ножницы десневые
 ■ Ciseaux a gencives
 ■ Zahnfleischscheren
 ■ Ножиці для ясен
 ■ Tijeras para encias
 ■ مقص صمغ



Landolt Operating Scissors

Deaver Operating,
Dissecting Scissors



Longe Operating
Dissecting Scissors

Incision Operating Scissors

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



Mixer Operating
Dissecting Scissors



TC28-247
15.5 cm
B/B Str.



TC28-248
15.5 cm
B/B Cvd.



Lexer Operating
Dissecting Scissors



TC28-249
16 cm
B/B Str.



TC28-251
16 cm
B/B Cvd.



Jorgenson Operating
Dissecting Scissors
TC28-253 22 cm
B/B, Full Cvd.



Richter Operating
Dissecting Scissors
TC28-254 14.5 cm
Angled



Fergusson Operating
Dissecting Scissors
TC28-255 14 cm
TC28-256 16 cm
TC28-257 18 cm
B/B, Angled

 Ножицы десневые
 Ciseaux à gencives
 Zahnfleischscheren
 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ



Wertheim Operating Dissecting Scissors
TC 28-258 14.5 cm
TC 28-259 22.5 cm
 B/B, Cvd.



Metzenbaum Dissecting Scissors
TC 28-260 20 cm
TC 28-261 23 cm
TC 28-262 28 cm
 B/B, Cvd.



Metzenbaum-Jamison Scissors
TC 28-263 15 cm
 S/S, Str.



TC 28-264
 15 cm
 B/B Str.

TC 28-265
 15 cm
 B/B Cvd.



Mayo Dissecting Scissors



TC 28-266
 14.5 cm
TC 28-267
 15 cm
TC 28-268
 17 cm
TC 28-269
 23 cm
 B/B, Str.



TC 28-270
 14.5 cm
TC 28-271
 15 cm
TC 28-272
 17 cm
TC 28-273
 23 cm
 B/B, Cvd.



Mayo-Stille Dissecting Scissors



TC 28-274
 15 cm
TC 28-275
 17 cm
TC 28-276
 19.5 cm
TC 28-277
 21 cm
 B/B, Str.



TC 28-278
 15 cm
TC 28-279
 17 cm
TC 28-280
 19.5 cm
TC 28-281
 21 cm
 B/B, Cvd.

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



TC28-282
22 cm
B/B, Str.

Mayo-Harrington Dissecting Scissors



TC28-283
22 cm
B/B, Cvd.



TC28-284
17 cm
B/B, Str.



TC28-285
17 cm
B/B, Cvd.

Mayo-Noble Dissecting Scissors



Metzenbaum Operating Dissecting Scissors
Blunt/Blunt Str.



TC28-286 14.5 cm
TC28-287 18 cm
TC28-288 20 cm
TC28-289 23 cm
TC28-290 25 cm
TC28-291 28 cm
B/B Str.



TC28-292 14.5 cm
TC28-293 18 cm
TC28-294 20 cm
TC28-295 23 cm
TC28-296 25 cm
TC28-297 28 cm
B/B Cvd.



Metzenbaum-Nelson Scissors



TC28-298 18 cm
TC28-299 20 cm
TC28-300 23 cm
TC28-301 25 cm
TC28-302 28 cm
TC28-303 30 cm
B/B Str.



TC28-304 18 cm
TC28-305 20 cm
TC28-306 23 cm
TC28-307 25 cm
TC28-308 28 cm
TC28-309 30 cm
S/B Str.



TC28-310 18 cm
TC28-311 20 cm
TC28-312 23 cm
TC28-313 25 cm
TC28-314 28 cm
TC28-315 30 cm
S/S Str.



Metzenbaum-Nelson Scissors



TC28-316 18 cm
TC28-317 20 cm
TC28-318 23 cm
TC28-319 25 cm
TC28-320 28 cm
TC28-321 30 cm
B/B Cvd.



TC28-322 18 cm
TC28-323 20 cm
TC28-324 23 cm
TC28-325 25 cm
TC28-326 28 cm
TC28-327 30 cm
S/B Cvd.



TC28-328 18 cm
TC28-329 20 cm
TC28-330 23 cm
TC28-331 25 cm
TC28-332 28 cm
TC28-333 30 cm
S/S Cvd.

 Ножницы десневые
 Ciseaux à gencives
 Zahnfleischscheren
 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ



Metzenbaum-Fino Scissors

TC28-334
14.5 cm
B/B Str.

TC28-335
14.5 cm
S/S Str.


TC28-336
14.5 cm
B/B Cvd.

TC28-337
14.5 cm
S/S Cvd.

TC28-338
18 cm
TC28-339
20 cm
TC28-340
23 cm
B/B Str.

TC28-341
18 cm
TC28-342
20 cm
TC28-343
23 cm
B/B Cvd.


TC28-344
18 cm
TC28-345
20 cm
TC28-346
23 cm
B/B Slightly Cvd.




Metzenbaum-Fino Scissors



Kaye Face Lift Scissors
TC28-347 11.5 cm
TC28-348 14 cm
TC28-349 15 cm
TC28-350 18 cm
TC28-351 20 cm
B/B, Str.



Kaye Face Lift Scissors
TC28-352 11.5 cm
TC28-353 14 cm
TC28-354 15 cm
TC28-355 18 cm
TC28-356 20 cm
B/B, Cvd.



Gorney Face Lift Scissors
TC28-357 18 cm
TC28-358 18 cm serrated
TC28-359 23 cm serrated
B/B, Str.




TC28-360 18 cm
TC28-361 23 cm
B/B, Cvd, Serrated



TC28-362 18 cm
B/B, Cvd.


Ножницы десневые
Ciseaux a gencives
Zahnfleischschere

Ножиці для ясен
Tijeras para encias
مقص صمغ




Rees Face Lift Scissors
TC28-363 17 cm
B/B, Cvd.

Rees Face Lift Scissors
TC28-364 17 cm
B/B, Str.



Kilner Operating Dissecting Scissors
TC28-365 12 cm
TC28-366 15 cm
TC28-367 18 cm
TC28-368 20 cm
B/B Str.

Kilner Operating Dissecting Scissors
TC28-369 12 cm
TC28-370 15 cm
TC28-371 18 cm
TC28-372 20 cm
B/B Cvd.



Reynolds Operating, Dissecting Scissors
TC28-373 15 cm
TC28-374 18 cm
S/S Str.



Reynolds Operating, Dissecting Scissors
TC28-375 15 cm
TC28-376 18 cm
S/S Cvd.



Sanvenero Operating Dissecting Scissors

28-377
14 cm
S/S Str.

28-378
14 cm
S/S Cvd.



Dubois Dissecting & Gynecological Scissors
TC28-379 27 cm
B/B Str.



Dubois Dissecting & Gynecological Scissors
TC28-380 27 cm
B/B Cvd.



Neumann Vascular Scissors
TC28-381 12.5 cm
S/S Double Cvd.

■ Ножницы десневые
 ■ Ciseaux a gengivas
 ■ Zahnfleischscheren
 ■ Ножницы для ясен
 ■ Tijeras para encias
 ■ مقص صمغ



Killner Tendon Scissors
TC28-382 11.5 cm
Straight, B/B

Littler Tendon Scissors

Strabismus Tendon Scissors

Potts-Reynolds Tenotomy Scissors
TC28-387 12 cm
S/S Cvd.

TC28-383 11.5 cm
 TC28-384 14 cm
 Cvd, S/S

TC28-385 11.5 cm
 Str, B/B
 TC28-386 11.5 cm
 Cvd, B/B



Wagner Gum Scissors
Plastic surgery

Brophy (sullivan) Scissors

Potts Tenotomy Scissors
TC28-392 15 cm
TC28-393 17 cm
TC28-394 20 cm
S/S

Wilmer Iris Scissors
TC28-395 9 cm
S/S Angled

TC28-388
 12 cm
 S/S Str.


TC28-390
 14.5 cm
 S/S Str.

TC28-389
 12 cm
 S/S Cvd.


TC28-391
 14.5 cm
 S/S Cvd.

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren


Ножиці для ясен
Tijeras para encias
مقص صمغ




**Potts-Smith
Vascular Scissors
S/S Angled**




**TC28-396
19 cm**



**TC28-397
19 cm**



**TC28-398
19 cm**



**TC28-399
19 cm**



**Potts-Smith Scissors
TC28-400 19.5 cm
B/B Angled.**



**Cardiovascular Scissors
TC28-401 12 cm
Probe End**



**Enterotomy Scissors
TC28-402 21 cm**



**Enterotomy Scissors
TC28-403 21 cm**



**Enterotomy Scissors
TC28-404 21 cm**






**Spencer Stitch Scissors
TC28-405 9 cm
TC28-406 11 cm
TC28-407 14 cm
Str.**

Ножницы десневые
Ciseaux à gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encías
مقص صمغ



Artery Scissors
Probe end Angled




Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren


Ножиці для ясен
Tijeras para encias
مقص صمغ




Martin-Cartilage Scissors
TC28-418 20 cm
S/S Serrated



Mc Indoe Scissors
TC28-419 19 cm
S/S Str.



Potts-Smith Scissors
TC28-420 19 cm
B/B Cvd.




Crafoord Scissors
TC28-421 30 cm
B/B Cvd.



Stevens Scissors

TC28-422
10 cm
TC28-423
11 cm
S/S Str.

TC28-424
10 cm
TC28-425
11 cm
S/S Cvd.



Stevens Scissors

TC28-426
10 cm
TC28-427
11 cm
B/B Str.

TC28-428
10 cm
TC28-429
11 cm
B/B Cvd.



Graefe Scissors
TC28-430 10.5 cm
B/B Angled

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



Cottle Scissors
TC28-431 16 cm
B/B, Angled





Fomon Scissors
TC28-432 14.5 cm
B/B, Angled



Sailer Scissors
TC28-433 17 cm
B/B, Angled



Harvey Scissors
TC28-434 12.5 cm
Saw Edge



Schoemaker-Loth Scissors
TC28-435 12.5 cm



Eiselsberg Scissors
TC28-436 11 cm

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren


Ножиці для ясен
Tijeras para encias
مقص صمغ




Scissors
TC28-437 9 cm
TC28-438 10.5 cm
S/S, Str.



Scissors
TC28-439 13 cm
TC28-440 15 cm
S/B, Str.



Scissors
TC28-441 9 cm
TC28-442 10.5 cm
S/S, Str.



Scissors
TC28-443 9 cm
TC28-444 10.5 cm
S/S, Cvd.




Cuticle Scissors



TC28-445
10 cm
TC28-446
11.5 cm
S/S, Str.



TC28-447
10 cm
TC28-448
11.5 cm
S/S, Cvd.



Cuticle Scissors



TC28-449
9 cm
TC28-450
10 cm
TC28-451
11.5 cm
S/S, Str.



TC28-452
9 cm
TC28-453
10 cm
TC28-454
11.5 cm
S/S, Cvd.

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقصص صمغ



Castroviejo Scissors
Spring Type




TC28-455
9 cm
S/B Str.



TC28-456
9 cm
S/B Cvd.




TC28-457
9 cm
S/S Str.



TC28-458
9 cm
S/S Cvd.




TC28-459
11.5 cm
B/B Str.



TC28-460
11.5 cm
B/B Cvd.




TC28-461
11.5 cm
S/B Str.



TC28-462
11.5 cm
S/B Cvd.




TC28-463
11.5 cm
S/S Str.



TC28-464
11.5 cm
S/S Cvd.




TC28-465
14 cm
B/B Str.



TC28-466
14 cm
B/B Cvd.




TC28-467
14 cm
S/B Str.



TC28-468
14 cm
S/B Cvd.



TC28-469
14 cm
S/S Str.



TC28-470
14 cm
S/S Cvd.

Ножницы десневые
Ciseaux a gengivas
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



Noyes Scissors
Spring Handle



Noyes Scissors
Spring Handle
TC28-474 12 cm
Angled,



Wecker Scissors
Angled



Micro Scissors
TC28-478 12 cm
S/S, Cvd.



Westcott Micro Scissors
TC28-479 11 cm
Cvd-Micro S/S



Castroviejo Corneal Micro Scissors
TC28-480 11.5 cm
Left Angled,
Spring Handle



Castroviejo Corneal Micro Scissors
TC28-481 11.5 cm
Right Angled,
Spring Handle



Troutman Micro Scissors
TC28-482 10 cm
Left Angled,
Spring Handle



Troutman Micro Scissors
TC28-483 10 cm
Right Angled,
Spring Handle


■ Ножницы десневые
 ■ Ciseaux a gencives
 ■ Zahnfleischscheren
 ■ Ножницы для ясен
 ■ Tijeras para encias
 ■ مقص صمغ



Westcott Micro Scissors




Westcott Micro Scissors



14 cm




Dissection Micro Scissors
TC 28-488 14 cm
B/B, Cvd.




Adventitia Micro Scissors
TC 28-489 14 cm
S/S, Str.




Delicate Micro Scissor




Castroviejo Micro Scissors



18 cm




Dissection Micro Scissors
TC 28-494 18 cm B/B Cvd.




Adventitia Micro Scissors
TC 28-495 18 cm S/S Str.

■ Ножницы десневые
 ■ Ciseaux a gencives
 ■ Zahnfleischscheren

■ Ножницы для ясен
 ■ Tijeras para encias
 ■ مقص صمغ



15.0 cm, 6"



Dissection Micro Scissors
TC28-496 15 cm
S/S, Cvd.



Adventitia Micro Scissors
TC28-497 15 cm
S/S, Str.



Dissection Micro Scissors
TC28-498 15 cm
S/S, Str.



Adventitia Micro Scissors
TC28-499 15 cm
S/S, Cvd.



Dissection Micro Scissors
TC28-500 15 cm
S/S Cvd.



Adventitia Micro Scissors
TC28-501 15 cm
S/S, Str.

-Nerve-



Dissecting Scissors



TC28-502
11 cm
S/S, Str.



TC28-503
11 cm
S/S, Cvd.



TC28-504
18.5 cm
S/S, Str.



TC28-505
18.5 cm
S/S, Cvd.




Micro Scissors

■ Ножницы десневые
 ■ Ciseaux a gengivas
 ■ Zahnfleischscheren

■ Ножиці для ясен
 ■ Tijeras para encias
 ■ مقصص صمغ



Yasargill Scissors



TC28-506 16 cm
TC28-507 20 cm
TC28-508 22 cm
 S/S, Str.



TC28-509 16 cm
TC28-510 20 cm
TC28-511 22 cm
 S/S, Cvd.



TC28-512 16 cm
TC28-513 20 cm
TC28-514 22 cm
 S/S, Str.



TC28-515 16 cm
TC28-516 20 cm
TC28-517 22 cm
 S/S, Cvd.



TC28-518 16 cm
TC28-519 20 cm
TC28-520 22 cm
 B/B, Str.



TC28-521 16 cm
TC28-522 20 cm
TC28-523 22 cm
 B/B, Cvd.



Micro Scissors
S/S, Angled




TC28-524
 18 cm



TC28-525
 18 cm



TC28-526
 18 cm



TC28-527
 18 cm



TC28-528
 18 cm




Micro Scissors
S/S, Angled pointed probe



TC28-529
 18 cm



TC28-530
 18 cm



TC28-531
 18 cm



TC28-532
 18 cm



TC28-533
 18 cm

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren


Ножиці для ясен
Tijeras para encias
مقص صمغ



Bruns Scissors

TC28-534
24 cm
Cvd, Probe End

TC 28-535
24 cm
Cvd, Serrated
Probe End



**Lister Bandage
Plaster Scissors**
TC28-536 9cm
TC28-537 11 cm
TC28-538 14 cm
TC28-539 16 cm
TC28-540 18 cm
TC28-541 20 cm
TC28-542 24 cm




**Burnham Bandage
Plaster Scissors**
TC28-543 14 cm
Angled



**Braun-Stadler Plaster
Scissors**
TC28-544 14 cm
TC28-545 20 cm
B/B, Angled



Busch Plaster Scissors
TC28-546 13 cm
TC28-547 16 cm
B/B, Cvd.



Modified-USA Plaster Scissors
TC 28-548 10 cm



Seutin Plaster Scissors
TC28-549 20 cm
Probe End



Seutin Plaster Scissors
TC28-550 23 cm
Probe End

■ Ножницы десневые
 ■ Ciseaux a gencives
 ■ Zahnfleischscheren
 ■ Ножиці для ясен
 ■ Tijeras para encias
 ■ مقص صمغ

T.C. / CERAMIC COATED SCISSORS



Mayo Dissecting Scissors

TC 28-551
 14.5cm
TC28-552
 17cm
TC28-553
 23cm
 B/B, Str.

TC28-554
 14.5cm
TC28-555
 17cm
TC 28-556
 23cm
 B/B, Cvd.



Metzenbaum Standard Scissors

TC28-557
 14.5cm
TC28-558
 18cm
TC28-559
 20cm
TC28-560
 23cm
TC28-561
 26cm
TC28-562
 28cm
TC28-563
 30cm
 B/B, Str.

TC 28-564
 14.5cm
TC 28-565
 18cm
TC 28-566
 20cm
TC 28-567
 23cm
TC 28-568
 26cm
TC 28-569
 28cm
TC 28-570
 30cm
 B/B, Cvd.



Locklin Dissecting Scissors, S/S

TC28-571
 16cm
 1 Blade
 Serrated



Locklin Dissecting Scissors, S/S

TC 28-572
 16cm
 1 Blade
 Serrated



Dean Dissecting Scissors, S/S


TC 28-573
 17cm
 1 Blade
 Serrated

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



TC28-574
Iris
11cm
Straight



TC28-575
Iris
11cm
Curved




TC28-576
Goldman-Fox
13cm
Saw edge, straight



TC28-577
Goldman-Fox
13cm
Saw edge, Curved



TC28-578
Iris Double Bend
11cm



TC28-579
La Grange
11.5cm
Saw edge,
double Curved



TC28-580
Dean Scissors TC
17cm curve



TC28-581
Kelly
16cm
Straight



TC28-582
Kelly
16cm
Curved

 Ножицы десневые
 Ciseaux à gencives
 Zahnfleischscheren

 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ




**Goldmann-Fox
Operating Scissors
T/C Insert**



TC28-583
13 cm
S/S Str.



TC28-584
13 cm
S/S Cvd.




**Lister Bandage Plaster
Scissors T/C Insert**
one large ring
TC 28-585 14 cm
TC 28-586 16 cm
TC 28-587 18 cm
TC 28-588 20 cm




**Lister Bandage Plaster
Scissors, T/C Insert**
TC 28-589 9 cm
TC 28-590 11.5 cm
TC 28-591 14 cm
TC 28-592 16 cm
TC 28-593 18 cm
TC 28-594 20 cm



TC28-595
Universal
12 cm



TC28-596
Beebee
11 cm
Saw edge,
straight



TC28-597
Beebee
11 cm
Saw edge,
curved



TC 28-598
Beebee
12 cm
Saw edge,
straight




TC28-599
Beebee
12 cm
Saw edge,
curved

Ножницы десневые
Ciseaux a gencives
Zahnfleischscheren

Ножиці для ясен
Tijeras para encias
مقص صمغ



TC28-600
Goldman Fox Scissors
12.5cm Curved



TC28-601
Iris Scissors
11.5cm Curved
Serrated For Slip-Free Cutting



TC28-602
LaGrange Scissors
11.5cm Double Curved
Serrated For Slip-Free Cutting



TC28-603
Castroviejo Micro Scissors
12cm Curved




TC28-604
Super Cut Goldman Fox
13.2cm Curved
Serrated For Slip-Free Cutting



TC28-605
Super Cut Iris
11cm Straight
Serrated For Slip-Free Cutting



TC28-606
Super Cut Iris
11cm Curved
Serrated For Slip-Free Cutting



TC28-607
Super Cut LaGrange
11.5cm
Double Curved
Serrated For Slip-Free Cutting


 Ножицы десневые
 Ciseaux a gencives
 Zahnfleischscheren

 Ножиці для ясен
 Tijeras para encias
 مقصص صمغ



Operating Scissors, S/S Str. T/C Insert

TC 28-608	11.5 cm
TC28-609	13 cm
TC28-610	14 cm
TC28-611	15 cm
TC28-612	16 cm
TC28-613	17.5 cm
TC28-614	18.5 cm
TC28-615	20 cm




Operating Scissors, S/S Cvd. T/C Insert

TC 28-616	11.5 cm
TC 28-617	13 cm
TC 28-618	14 cm
TC 28-619	15 cm
TC 28-620	16 cm
TC 28-621	17.5 cm
TC 28-622	18.5 cm
TC 28-623	20 cm



Operating Scissors, B/B Str. T/C Insert

TC 28-624	11.5 cm
TC28-625	13 cm
TC28-626	14 cm
TC28-627	15 cm
TC28-628	16 cm
TC28-629	17.5 cm
TC28-630	18.5 cm
TC28-631	20 cm



Operating Scissors, B/B Cvd. T/C Insert

TC 28-632	11.5 cm
TC28-633	13 cm
TC28-634	14 cm
TC28-635	15 cm
TC28-636	16 cm
TC28-637	17.5 cm
TC28-638	18.5 cm
TC28-639	20 cm



Operating Scissors, S/B Standard Str. T/C Insert

TC 28-640	11.5 cm
TC28-641	13 cm
TC28-642	14 cm
TC28-643	15 cm
TC28-644	16 cm
TC28-645	17.5 cm
TC28-646	18.5 cm
TC28-647	20 cm



Operating Scissors, S/B Standard Cvd. T/C Insert

TC 28-648	11.5 cm
TC 28-649	13 cm
TC 28-650	14 cm
TC 28-651	15 cm
TC 28-652	16 cm
TC 28-653	17.5 cm
TC 28-654	18.5 cm
TC 28-655	20 cm

 Ножицы десневые
 Ciseaux à gencives
 Zahnfleischscheren

 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ



Metzenbaum Operating Scissors T/C Insert

TC28-656
 11.5 cm
TC28-657
 14.5 cm
 B/B Str.

TC 28-658
 11.5 cm
TC 28-659
 14.5 cm
 S/B Str.

TC28-660
 11.5 cm
TC28-661
 14.5 cm
 S/S Str.

TC28-662
 11.5 cm
TC28-663
 14.5 cm
 B/B Cvd.

TC28-664
 11.5 cm
TC28-665
 14.5 cm
 S/B Cvd.

TC 28-666
 11.5 cm
TC 28-667
 14.5 cm
 S/S Cvd.



Metzenbaum Operating Scissors T/C Insert



Metzenbaum-Nelson Operating Scissors T/C Insert

TC 28-668
 18 cm
TC 28-669
 20 cm
TC 28-670
 23 cm
TC 28-671
 25 cm
TC 28-672
 28 cm
TC 28-673
 30 cm
 B/B Str.

TC28-674
 18 cm
TC28-675
 20 cm
TC28-676
 23 cm
TC28-677
 25 cm
TC28-678
 28 cm
TC28-679
 30 cm
 S/B Str.

TC28-680
 18 cm
TC28-681
 20 cm
TC28-682
 23 cm
TC28-683
 25 cm
TC28-684
 28 cm
TC28-685
 30 cm
 S/S Str.

TC 28-686
 18 cm
TC 28-687
 20 cm
TC 28-688
 23 cm
TC 28-689
 25 cm
TC 28-690
 28 cm
TC 28-691
 30 cm
 B/B Cvd.

TC 28-692
 18 cm
TC 28-693
 20 cm
TC 28-694
 23 cm
TC 28-695
 25 cm
TC 28-696
 28 cm
TC 28-697
 30 cm
 S/B Cvd.

TC 28-698
 18 cm
TC 28-699
 20 cm
TC 28-700
 23 cm
TC 28-701
 25 cm
TC 28-702
 28 cm
TC 28-703
 30 cm
 S/S Cvd.



Metzenbaum-Nelson Operating Scissors T/C Insert

 Ножицы десневые
 Ciseaux à gencives
 Zahnfleischscheren

 Ножиці для ясен
 Tijeras para encías
 مقصص صمغ



Metzenbaum-Fino
Operating Scissors
T/C Insert

TC28-704 14.5 cm
TC28-705 18 cm
TC28-706 20 cm
TC28-707 23 cm
 B/B Str.



TC28-708 14.5 cm
TC28-709 18 cm
TC28-710 20 cm
TC28-711 23 cm
 B/B Cvd.



TC28-712 14.5 cm
TC28-713 18 cm
TC28-714 20 cm
TC28-715 23 cm
 S/S Str.



TC28-716 14.5 cm
TC28-717 18 cm
TC28-718 20 cm
TC28-719 23 cm
 S/S Cvd.



Metzenbaum-Fino
Operating Scissors
T/C Insert



Mayo Dissecting &
Gynecological Scissors
T/C Insert

TC28-720
14.5 cm
TC28-721
15 cm
TC28-722
17 cm
TC28-723
23 cm
 B/B, Str.



TC28-724
14.5 cm
TC28-725
15 cm
TC28-726
17 cm
TC28-727
23 cm
 B/B, Cvd.



TC28-728
15 cm
TC28-729
17 cm
TC28-730
19.5 cm
TC28-731
21 cm
 B/B, Str.




TC28-732
15 cm
TC28-733
17 cm
TC28-734
19.5 cm
TC28-735
21 cm
 B/B, Cvd.



Mayo-stille Dissecting &
Gynecological Scissors
T/C Insert


■ Ножницы десневые
 ■ Ciseaux a gencives
 ■ Zahnfleischscheren
 ■ Ножиці для ясен
 ■ Tijeras para encias
 ■ مقص صمغ



TC 28-736
14 cm
B/B Str.

TC 28-737
14 cm
B/B Cvd.

TC 28-738
14 cm
S/B Str.



TC 28-739
14 cm
S/B Cvd.



TC 28-740
14 cm
S/S Str.



TC 28-741
14 cm
S/S Cvd.



TC 28-742
14 cm
S/S Str.



TC 28-743
14 cm
S/S Cvd.



Joseph Operating
Scissors T/C Insert




Reynolds Operating
Scissors T/C Insert



TC 28-744
16 cm
TC 28-745
18 cm
S/S Str.



TC 28-746
16 cm
TC 28-747
18 cm
S/S Cvd.



Foman Operating
Scissors T/C Insert
TC 28-748
13 cm Full Bent



Stevens Operating
Scissors T/C Insert



TC 28-749
10 cm
TC 28-750
11 cm
Str.




TC 28-751
10 cm
TC 28-752
11 cm
Cvd.

Зажимы артериальные
Pinces hemostatiques
Arterienklemmen

Затискачі артеріальні
Pinzas hemostaticas
ملقطه هایموستاتیک



Hartmann (Baby Mosquito)
TC28-753 10 cm straight
TC28-754 10 cm curved



Halstead Mosquito
TC28-755 10cm
TC28-756 12.5cm
TC28-757 14cm
straight
TC28-758 10cm
TC28-759 12.5cm
TC28-760 14cm
curved



Halstead Mosquito 1x2
TC28-761 10cm
TC28-762 12.5cm
TC28-763 14cm
straight
TC28-764 10cm
TC28-765 12.5cm
TC28-766 14cm
Curved



Kelly
TC28-767 14cm
TC28-768 16cm
TC28-769 18cm
TC28-770 20cm
TC28-771 22cm
TC28-772 24cm
straight
TC28-773 14cm
TC28-774 16cm
TC28-775 18cm
TC28-776 20cm
TC28-777 22cm
TC28-778 24cm
curved



Rochester-Pean
TC28-779 14cm
TC28-780 16cm
TC28-781 18cm
TC28-782 20cm
TC28-783 22cm
TC28-784 24cm
TC28-785 26cm
TC28-786 28cm
TC28-787 30cm
straight
TC28-788 14cm
TC28-789 16cm
TC28-790 18cm
TC28-791 20cm
TC28-792 22cm
TC28-793 24cm
TC28-794 26cm
TC28-795 28cm
TC28-796 30cm
curved





Rochester-Kocher 1x2
TC28-797 14cm
TC28-798 16cm
TC28-799 18cm
TC28-800 20cm
TC28-801 22cm
TC28-802 24cm
TC28-803 26cm
TC28-804 28cm
TC28-805 30cm
straight
TC28-806 14cm
TC28-807 16cm
TC28-808 18cm
TC28-809 20cm
TC28-810 22cm
TC28-811 24cm
TC28-812 26cm
TC28-813 28cm
TC28-814 30cm
curved



Kelly-Rankin
TC28-815 16 cm
straight
TC28-816 16 cm
curved

■ Зажимы артериальные ■ Затискачі артеріальні
 ■ Pinces hemostatiques ■ Pinzas hemostáticas
 ■ Arterienklemmen ■ ملقط هایموستاتیک



Зажимы артериальные
Pinces hemostatiques
Arterienklemmen


Затискачі артеріальні
Pinzas hemostaticas
ملقط هایموستاتیک




Pean Haemostatic Forceps
Slender Type




TC 28-854 13 cm
TC 28-855 14 cm
straight




TC 28-856 13 cm
TC 28-857 14 cm
curved



Spencer-Wells Haemostatic Forceps
Slender Type




TC 28-858 13 cm
TC 28-859 15 cm
TC 28-860 17.5 cm
TC 28-861 20 cm
straight




TC 28-862 13 cm
TC 28-863 15 cm
TC 28-864 17.5 cm
TC 28-865 20 cm
curved



Kocher Haemostatic
Forceps Strong Pattern
1x2 teeth



TC 28-866 13cm
TC 28-867 14cm
Straight



TC 28-868 13cm
TC 28-869 14cm
curved



Rochester-Carmalt
Haemostatic Forceps



TC 28-870 16 cm
TC 28-871 20 cm
Straight



TC 28-872 16 cm
TC 28-873 20 cm
Curved

■ Зажимы артериальные ■ Затискачі артеріальні
 ■ Pinces hemostatiques ■ Pinzas hemostáticas
 ■ Arterienklemmen ■ ملقط هایموستاتیک



**Jones
Towel Forceps**
TC 28-874 9 cm
Slight Curve



**Engl. Modell
Towel Forceps**
TC 28-875 9 cm




**Schaedel Towel
Forceps**
TC 28-876 9 cm



Mayo
TC 28-877 14cm



Backhaus Clip Towel Forceps
TC 28-878 11 cm
with clip for tubing or cable



Roeder Towel Forceps
TC 28-879 13 cm
TC 28-880 15 cm
Curved



Bernhard Towel Forceps
TC 28-881 16.5 cm

Зажимы артериальные
Pinces hemostatiques
Arterienklemmen

Затискачі артеріальні
Pinzas hemostaticas
ملقط هایموستاتیک




Pratt Haemostatic Forceps
TC28-882 15 cm



Pennington / Siezing Haemostatic Forceps
TC28-883 15.5 cm
TC28-884 20 cm
Straight




TC28-885 12.5 cm
TC28-886 14 cm
TC28-887 16 cm



Pean Forceps


TC28-888 16 cm
TC28-889 18 cm



Lister Sinus Forceps

TC28-890 13cm
TC28-891 15.5cm
TC28-892 17.5cm

■ Зажимы артериальные ■ Затискачі артеріальні
 ■ Pinces hemostatiques ■ Pinzas hemostáticas
 ■ Arterienklemmen ■ ملقط هایموستاتیک



Foerster-Ballenger Sponge Forceps
& (bullet holding forceps)




TC28-893 18 cm
TC28-894 20 cm
TC28-895 25 cm
 Straight Smooth



TC28-896 18 cm
TC28-897 20 cm
TC28-898 25 cm
 Straight Serrated



TC28-899 18 cm
TC28-900 20 cm
TC28-901 25 cm
 Curved Smooth



TC28-902 18 cm
TC28-903 20 cm
TC28-904 25 cm
 Curved Serrated

STERILIZING FORCEPS



TC28-905 20 cm
TC28-906 25 cm
TC28-907 30 cm



Rogge
TC28-908 20 cm
TC28-909 30 cm




Cheatel
 Box Joint
TC28-910
 27 cm




Cheatel
 Screw Joint
TC28-911
 27 cm

Зажимы артериальные
Pinces hemostatiques
Arterienklemmen

Затискачі артеріальні
Pinzas hemostáticas
ملقط هائموستاتیک




Backhaus
TC28-912
Not-perforating and
not-damaging 14cm



Backhaus
TC28-913 11cm
TC28-914 13cm
TC28-915 14cm
TC28-916 15cm



Backhaus
TC28-917
9cm



Gross Maier Dressing Forceps
TC28-918 16cm straight
TC28-919 20cm straight
TC28-920 25cm straight

TC28-921
16cm
TC28-922
20cm
TC28-923
25cm
curved



TC28-924
Young
Tongue Forceps
17cm





TC28-925
Collins
Tongue Forceps
17cm

Spare rubber


Стерилизационные и тканевые зажимы
Stérilisation et pinces pour tissus
Esterilizadores y sujetadores de tejidos

Стерилізаційні та тканинні затискачі
Esterilização e grampos de tecido
تعقيم و المشبك الأنسجة



Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачі
Porta-agujas
اصحاب ابرة




Castroviejo


TC28-937
13 cm

TC28-938
15 cm

TC28-939
18 cm




TC28-940
Halsey
13 cm



Baumgartner

TC28-941 13cm
TC28-942 14cm
TC28-943 16cm
TC28-944 18cm
TC28-945 20cm
TC28-946 25cm
TC28-947 30cm




Derf

TC28-948
12 cm


TC28-949
12 cm




TC28-950
Delicate
15cm



TC28-951
Adson
18cm



Mayo-Hegar
TC28-952 13cm
TC28-953 14cm
TC28-954 15cm
TC28-955 16cm
TC28-956 18cm
TC28-957 20cm
TC28-958 23cm
TC28-959 24cm
TC28-960 26cm
TC28-961 30cm
Serrated Grooved



Crile-Wood
TC28-962 15cm
TC28-963 18cm
TC28-964 20cm
TC28-965 22cm
TC28-966 25cm

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة




Webster
Needle Holder

Heaney Needle Holders

Crile-Murray
Needle Holders

Doyen Needle Holders




Olsen-Hegar Needle
Holder With Scissors

Olsen-Hegar Needle
Holder With Scissors

TC28-996
Gillies Needle
Holder With Scissors
16cm

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачі
Porta-agujas
اصحاب ابرة



Kilner Needle Holders

TC28-997
13 cm
serrated



Crile
TC28-998
15 cm



Toennis
TC28-999
18 cm



Mathieu

TC28-1000 14 cm
TC28-1001 17 cm

TC28-1002
20 cm



Mathieu

TC28-1003 14 cm
TC28-1004 17 cm

TC28-1005
20 cm



Mathieu

TC28-1006 14 cm
TC28-1007 17 cm

TC28-1008
20 cm



Mathieu

TC28-1009 14 cm
TC28-1010 17 cm

TC28-1011
20 cm

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة



TC28-1012
14 cm
Curved, Smooth
without lock
curved



TC28-1013
14 cm
Straight, Smooth
without lock
straight



TC28-1014
14 cm
Curved
Smooth
with lock

TC28-1015
14 cm
Curved
Smooth
without lock

TC28-1016
14 cm
Straight
Smooth
with lock

TC28-1017
14 cm
Straight
Smooth
without lock

TC28-1018
18 cm
Curved
Smooth
with lock

TC28-1019
18 cm
Curved
Smooth
without lock

TC28-1020
18 cm
Straight
Smooth
with lock

TC28-1021
18 cm
Straight
Smooth
without lock

Barraquer Needle Holders

Castroviejo Needle Holders



Castroviejo Olsen Heager /
Needle Holders / Scissors
TC28-1022 14 cm
TC28-1023 15 cm



TC28-1024
13 cm
with lock
straight

TC28-1025
13 cm
without lock
straight



TC28-1026
13 cm
with lock
curved

TC28-1027
13 cm
without lock
curved



TC28-1028
14 cm
Curved, Smooth
with lock curved

TC28-1029
14 cm
Straight, Smooth
with lock straight

Castroviejo Needle Holders

Barraquer Needle Holders

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة



TC28-1030
16 cm
TC28-1031
20 cm
TC28-1032
23 cm
Smooth Str.



TC28-1033
16 cm
TC28-1034
20 cm
TC28-1035
23 cm
Smooth Cvd.



TC28-1036
16 cm
Smooth Cvd.



TC28-1037
12 cm
Smooth Str
with lock



TC28-1038
12 cm
Smooth Cvd
with lock



TC28-1039
10.5 cm
Smooth
Str.



TC28-1040
10.5 cm
Smooth
Cvd.



TC28-1041
11.5 cm
Smooth
Str.



TC28-1042
11.5 cm
Smooth
Cvd.

Micro Needle Holders

Micro Needle Holders

Castroviejo Needle
Holders

Castroviejo Needle
Holders

Barraquer-Troutman
Needle Holders



TC28-1043
13.5 cm
Serrated



TC28-1044
13.5 cm
Serrated
Str.



TC28-1045
13.5 cm
Serrated
Cvd.



TC28-1046
14 cm



TC28-1047
17 cm

Mathieu-Kocher



TC28-1048
Boynton
12 cm

Kalt Needle Holders

Arruga Needle Holders

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачі
Porta-agujas
اصحاب ابرة


T.C GOLD




Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачі
Porta-agujas
اصحاب ابرة

T.C GOLD





Mathieu-Kocher




Mathieu



Olsen-Hegar Needle Holder With Scissors T.C



Olsen-Hegar Needle Holder With Scissors



TC28-1098
Gillies Needle Holder With Scissors
16cm



Castroviejo
TC28-1099 14cm straight
TC28-1100 18cm straight
TC28-1101 14cm curved
TC28-1102 18cm curved



Castroviejo
TC28-1103 14cm straight
TC28-1104 18cm straight
TC28-1105 14cm curved
TC28-1106 18cm curved

Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة



TC28-1107
Lichtenberg-
Ryder T.C
14 cm

TC28-1108
Lichtenberg-
Ryder T.C
14 cm

TC28-1109
Lichtenberg-
Ryder T.C
17 cm

TC28-1110
Lichtenberg-
Ryder T.C
17 cm



TC28-1111
Barraquer
18 cm

TC28-1112
Barraquer
18 cm

TC28-1113
Barraquer T.C
18 cm

TC28-1114
Barraquer T.C
18 cm



TC28-1115
Tönnis
18 cm



TC28-1116
Castroviejo
13 cm
Straight



TC28-1117
Kalt Needle Holders T.C
14 cm Straight




TC28-1118
Peet Splinter Forceps
11.5 cm



TC28-1119
La Grange
Curved

TC28-1120
All-Purpose Pliers

Иглодержатели
 Porte-aiguilles
 Nadelhalter
 Голкотримачі
 Porta-agujas
 اصحاب ابرة




Kilner
Needle Holders T.C

Wertheim Needle Holders T.C
TC28-1122 20 cm
TC28-1123 24 cm

Halsey Needle Holders T.C

Baumgartner Needle Holders T.C



Castroviejo Needle Holders T.C

Castroviejo Needle Holders T.C

Mathieu-Ryder Needle
Holders T.C

Crile

Инструменты для системы коффердам


Pinces porte-crampons a digue

Kofferdam Klammerzangen

Инструменти для системи коффердам

Pinzas de engrapar para dique de caucho

السند المطاط الملقط المشبك



TC28-1137
Brewer



TC28-1138
Stockes



TC28-1139
Palmer



TC28-1140
Ivory
light weight



TC28-1141
Anisworth




TC28-1142
Ivory




TC28-1143
Elastic
Seaparing & Ligature
Plier



TC28-1144
Elastic
Seaparing & Ligature
Plier



TC28-1145
Small
Rubber Dam Frame



TC28-1146
Medium
Rubber Dam Frame



TC28-1147
Large
Rubber Dam Frame

Резиновая дамба зажимы комплекты
Trousses de caoutchouc
Gummidammklammern

Гумова дамба затискачі комплекти
Conjuntos de abrazaderas de retención de caucho
المطاط البند المشابك مجموعات



Rubber Dam Kits




TC28-1148

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Rubber Dam Frame L
- 05 Rubber Dam Frame M
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7




TC28-1149

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Ivory
- 02 Ivory light weight
- 03 Clamps Stand
- 04 Rubber Dam Frame M
- 05 Rubber Dam Frame L
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 209
- 09 Rubber Dam Clamp Fig. 212
- 10 Rubber Dam Clamp Fig. 211
- 11 Rubber Dam Clamp Fig. 8
- 12 Rubber Dam Clamp Fig. 8A
- 13 Rubber Dam Clamp Fig. 14
- 14 Rubber Dam Clamp Fig. 14A
- 15 Rubber Dam Clamp Fig. 0
- 16 Rubber Dam Clamp Fig. 00
- 17 Rubber Dam Clamp Fig. 2
- 18 Rubber Dam Clamp Fig. 2A
- 19 Rubber Dam Clamp Fig. 7



TC28-1150
Rubber Dum Template

Резиновая дамба зажимы комплекты
Trousses de caoutchouc
Gummidammklammern

Гумова дамба затискачі комплекти
Conjuntos de abrazaderas de retención de caucho
المطاط البند المشبك مجموعات



TC28-1151

Stainless Steel Gold Plated Rubber Dam Clamp Tray with 11 Different Clamps
Fig. 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



TC 28-1152

Stainless Steel Rubber Dam Clamp Tray with 12 Different Clamps
Fig. 209, 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



TC28-1153

Stainless Steel Rubber Dam Clamp Tray with 8 Different Clamps
Fig. 1, 26, 29, 51, 60, 202, 205, 210



TC 28-1154

Stainless Steel Rubber Dam Clamp Tray
(without clamps)



TC28-1155
Surgical Kit



TC28-1156
Surgical Kit

Кламмера для коффердама
Crampons pour dique
Kofferdam Klammern

Кламмера для коффердама
Grapas para dique de caucho
المشايك السد المطاط



TC28-1157
Fig. 15



TC28-1158
Fig. 18



TC28-1159
Fig. 24



TC28-1160
Fig. 25



TC28-1161
Fig. 26



TC28-1162
Fig. 28



TC28-1163
Fig. 30



TC28-1164
Fig. 31



TC28-1165
Fig. 51



TC28-1166
Fig. 138



TC28-1167
Fig. 139



TC28-1168
Fig. 200



TC28-1169
Fig. 201



TC28-1170
Fig. 202



TC28-1171
Fig. 203



TC28-1172
Fig. 204



TC28-1173
Fig. 205



TC28-1174
Fig. 3



TC28-1175
Fig. W3



TC28-1176
Fig. 4



TC28-1177
Fig. 56A



TC28-1178
Fig. 7



TC28-1179
Fig. 7A



TC28-1180
Fig. 8



TC28-1181
Fig. 8A



TC28-1182
Fig. W8A



TC28-1183
Fig. 10



TC28-1184
Fig. 11



TC28-1185
Fig. 12A



TC28-1186
Fig. 13A



TC28-1187
Fig. 14



TC28-1188
Fig. 14A



TC28-1189
Fig. 14AN



TC28-1190
Fig. 22



TC28-1191
Fig. 27



TC28-1192
Fig. 207



TC28-1193
Fig. 29



TC28-1194
Fig. 206



TC28-1195
Fig. 208



TC28-1196
Fig. 209



TC28-1197
Fig. 1



TC28-1198
Fig. 1A



TC28-1199
Fig. 2



TC28-1200
Fig. 2A



TC28-1201
Fig. W2A



TC28-1202
Fig. 0



TC28-1203
Fig. 00



TC28-1204
Fig. 210



TC28-1205
Fig. 211



TC28-1206
Fig. 212



TC28-1207
Fig. 6



TC28-1208
Fig. 9



TC28-1209
Fig. N06



TC28-1210
Fig. U67



TC28-1211
Fig. P1



TC28-1212
Fig. P2



TC28-1213
Fig. G-1



TC28-1214
Fig. G-2

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменти для видалення коронок
Arranca-corona
مزيلات التاج

Complete With 3 Attachment



TC28-1215
Wynman Crown Gripper
No Lock T.C 16 cm



TC28-1216



TC28-1217
Crown Spreader, Crown Remover
Special Instrument for crown spreading
after separating, simple and safe
operation.



TC28-1218
Crown Remover



TC28-1219
Fig. 7



TC28-1220
Fig. 8



TC28-1221
Fig. 9



TC28-1222
Fig. 10



TC28-1223
Fig. 11



TC28-1224
Fig. 12



TC28-1225
Fig. 13
Tofflemire

Инструменты для удаления коронок

Arrache-couronnes

Teleskopkronenzangen


Инструменти для видалення коронок

Arranca-corona

مزيلات التاج



TC28-1226
Crown Spreaders



TC28-1227
Crown Spreaders
with Spring



TC28-1228
Crown Spreader



TC28-1229
Fig. 1



TC28-1230
Fig. 2



TC28-1231
Crown Spreaders
Fig. 3



TC28-1232
Crown Remover



TC28-1233
Crown Remover



TC28-1234
Crown Remover



Прочие товары
Divers produits et services
Sonstige Produkte

Інші товари
Varios Productos
المنتجات المتنوعة



TC 28-1235
X Ray Film Hanger



TC 28-1236
X Ray Film Hanger



TC 28-1237
X Ray Film Hanger



TC 28-1238
X Ray Film Hanger



TC 28-1239
Napkin Holder with
Plastic Chain



TC 28-1240
Napkin holder with metal chain



TC 28-1241
Mouth Props
Small



TC 28-1242
Mouth Props
Medium



TC 28-1243
Mouth Props
Large

Key Rings



TC 28-1244



TC 28-1245



TC 28-1246



TC 28-1247



TC 28-1248



TC 28-1249

Матрицедержатели
 Porte Matrices
 Matrizenspanner

Матрицеутримувачи
 Porta Matrices
 مصفوفة الخدم



TC28-1250
Fig. 5
Ivory 8
5 mm



TC28-1251
Fig. 6
Ivory 8
6 mm



TC28-1252
Fig. 7
Ivory 8
7 mm



TC28-1253
Ivory 8 N



TC28-1254
Didner 8 P



TC28-1255
Fig. 1
Tofflemire-Universal



TC28-1256
Fig. 2
Tofflemire-Junior



TC28-1257
Fig. 3
Tofflemire-Senior

Amalgam Well



TC28-1258
Stainless Steel
Amalgam Well
Medium



TC28-1259
Stainless Steel
Amalgam Well
large

■ Матрицедержатели ■ Матрицеутримувачи
 ■ Porte Matrices ■ Porta Matrices
 ■ Matrizenspanner ■ مصفوفة الخدم



TC28-1260
Tofflemire-Junior
45 mm



TC28-1261
Tofflemire-Universal
60 mm



TC28-1262
Matrix Ivory 8
6 mm



TC28-1263
Didner 8 P
70 mm



TC28-1264
Nyström
Left



TC28-1265
Nyström
Right



TC28-1266
Narrow



TC28-1267
Wide



TC28-1268
Ivory
with shields



TC28-1269
Ivory
without shields



TC28-1270
Ivory



TC28-1271
Elliot

■ Инструменты для амальгамы
 ■ Instruments pour amalgame
 ■ Amalgamis/strumente

■ Инструменты для амальгами
 ■ Instrumentos para amalgama
 ■ الصكوك الملتفم

Ø 2mm



TC28-1272
Fig. 10
16.5cm
English pattern

Ø 2mm



TC28-1273
Fig. 1
18.5cm
German pattern

Ø 2mm



TC28-1274
Fig. 1
15cm
Swiss pattern

Ø 2mm



TC28-1275
Fig. 1
15cm

Ø 2mm



TC28-1276
Fig. 2
15cm. 45°

Ø 2mm



TC28-1277
Fig. 3
14cm. 90°



TC28-1278
Chip Blower



TC28-1279



TC28-1280



TC28-1281



TC28-1282
Water Syringe

■ Инструменты для амальгамы ■ Инструменти для амальгами
 ■ Instruments pour amalgame ■ Instrumentos para amalgama
 ■ Amalgamis/strumente ■ الصكوك للملغم

AMALGAM CARRIERS

TIP SIZES FOR DISTAL OR SINGLE-END

1.5 mm

2.0 mm

2.8 mm

1.5 mm

2.0 mm

TC28-1283
Mini/
Single-End

TC28-1284
Regular
Single-End

TC28-1285
Large
Single-End

TC28-1286
Mini
Distal

TC28-1287
Regular
Distal

TIP SIZES FOR DOUBLE-END

1.5 mm

2.8 mm

3.2 mm

2.0 mm
TC28-1288
Mini/
Regular

2.0 mm
TC28-1289
Large/
Regular

2.0 mm
TC28-1290
Jumbo/
Regular

TC28-1291
Amalgam
Well

Single-End

TIP SIZES FOR DOUBLE-END

TC28-1292 - 1.5mm-2.0mm **TC28-1295** - 1.0mm-2mm
TC28-1293 - 2.0mm-2.8mm **TC28-1296** - 2.0mm-3mm
TC28-1294 - 2.0mm-3.2mm **TC28-1297** - 2.0mm-4mm


TC28-1298 - 1.5mm-2.0mm **TC28-1301** - 1.0mm-2mm
TC28-1299 - 2.0mm-2.8mm **TC28-1302** - 2.0mm-3mm
TC28-1300 - 2.0mm-3.2mm **TC28-1303** - 2.0mm-4mm

TC28-1304 - 1.5mm-2.0mm **TC28-1307** - 1.0mm-2mm
TC28-1305 - 2.0mm-2.8mm **TC28-1308** - 2.0mm-3mm
TC28-1306 - 2.0mm-3.2mm **TC28-1309** - 2.0mm-4mm

TC28-1310 - 1.5mm-2.0mm **TC28-1313** - 1.0mm-2mm
TC28-1311 - 2.0mm-2.8mm **TC28-1314** - 2.0mm-3mm
TC28-1312 - 2.0mm-3.2mm **TC28-1315** - 2.0mm-4mm

TC28-1316 - 1.5mm
TC28-1317 - 2.0mm
TC28-1318 - 2.0mm
TC28-1319 - 3.0mm

ROOT PLANNING CHISELES & BONE FILES



BONE FILES





Bone Splitting Chisels



Инструменты пародонтологические
Instruments pour la parodontologie
Instrumente für die Parodontologie

Інструменти пародонтологічні
Instrumentos para periodoncia
الصكوت اللثة



TC28-1349
1
Ochsenbin



TC28-1350
2
Ochsenbin



TC28-1351
Wakefield 1W
17cm



TC28-1352
Wakefield 2W
17cm



TC28-1353
Mead
19cm



TC28-1354
Mead
17cm



TC28-1355
Partson
18cm

ø 28 mm
Teflon
interchangeable
inserts
Autoclavable.

ø 26 mm
With
interchangeable
inserts
320 gr – 11 ¼ oz


ø 22 mm
Lead filled
200 gr – 7 oz




TC28-1356
Vickers
19 cm



TC28-1357
Herbert
19 cm




Parts 14 cm
TC28-1358 3mm
TC28-1359 4mm
TC28-1360 5mm
TC28-1361 6mm



Parts 17 cm
TC28-1362 3mm
TC28-1363 4mm
TC28-1364 5mm
TC28-1365 6mm



Parts 14 cm
TC28-1366 3mm
TC28-1367 4mm
TC28-1368 5mm
TC28-1369 6mm



Parts 17 cm
TC28-1370 3mm
TC28-1371 4mm
TC28-1372 5mm
TC28-1373 6mm

Инструменты пародонтологические

Instruments pour la parodontologie

Instrumente für die Parodontologie

Інструменти пародонтологічні

Instrumentos para periodoncia

الصكوت اللثة

PERIOSTEEL ELEVATORS



TC28-1374
Buser



TC28-1375
Ohl 10



TC28-1376
Howard 9H



TC28-1377
Allen 9A



TC28-1378
Prichard 3



TC28-1379
Seldin 23



TC28-1380



TC28-1381
Raspator, Raspatory
Length 17.5cm



TC28-1382



TC28-1383




TC28-1384
Raspator, Williger
Length 17cm, different widths
(4mm, 5mm or 7mm)



TC28-1385



TC28-1386
Trelat
14cm



Инструменты пародонтологические

Instruments pour la parodontologie

Instrumente für die Parodontologie

Інструменти пародонтологічні

Instrumentos para periodoncia

الصكوت اللثة



TC28-1387
Freer



TC28-1388
Freer 15



TC28-1389
Freer 16



TC28-1390
Goldman-Fox
14



TC28-1391
P 24G

6,0 mm



6,0 mm

TC28-1392
Molt
17.5cm

6,0 mm



4,5 mm

TC28-1393
Molt
17.5cm

6,5 mm



7,0 mm

TC28-1394
Meat
17.5cm

3,5 mm



3,5 mm

TC28-1395
Meat
17.5cm

2,5 mm



2,5 mm


TC28-1396
Meat
17.5cm

Шприцы стоматологические карпульные
Seringues à usage dentaire
Dentalspritzen

Шприці стоматологічні карпульні
Jeringas dentales
المحاقن الأسنان



Type A
TC28-1397 1.8ml
TC28-1398 2.2ml



Euro
Breech loading
TC28-1409 1.8ml
TC28-1410 2.2ml



Type C
TC28-1399 1.8ml
TC28-1400 2.2ml



Euro
Breech loading
TC28-1411 1.8ml
TC28-1412 2.2ml



Type C
TC28-1401 1.8ml
TC28-1402 2.2ml



Type Euro
TC28-1403 1.8ml
TC28-1404 2.2ml



Manual Aspirating Syringe
TC28-1413 1.8ml
TC28-1414 2.2ml



3-Ring
TC28-1405 1.8ml
TC28-1406 2.2ml



Quatro
TC28-1415 1.8ml
TC28-1416 2.2ml



2-Ring
TC28-1407 1.8ml
TC28-1408 2.2ml



Trio
TC28-1417 1.8ml
TC28-1418 2.2ml

Шприцы стоматологические карпульные
Seringues à usage dentaire
Dentalspritzen

Шприці стоматологічні карпульні
Jeringas dentales
المحاقن الأسنان



TC28-1419 1.8ml.
TC28-1420 2.2ml.
Pistol-Ject
Injects 0.20ml with every click
Lever action for high injection pressure

Autoclavable



TC28-1421 1.8ml
TC28-1422 2.2ml
Fig. 1001
Intraligamental
Anesthesia Syringe in
Black Colour Auto
Cleavable

Autoclavable



Ject Ultima
TC28-1423 1.8ml
TC28-1424 2.2ml
Injects 0,05 ml with every click
Supplied with angled and straight tubes



Swedish Type
TC28-1425 1.8ml
TC28-1426 2.2ml



Ject angled
TC28-1427 1.8ml
TC28-1428 2.2ml



Ject straight
TC28-1429 1.8ml
TC28-1430 2.2ml



Disc System
TC28-1431 1.8ml
TC28-1432 2.2ml



Disc System
TC28-1433 1.8ml
TC28-1434 2.2ml



Disc System
TC28-1435 1.8ml
TC28-1436 2.2ml



Projection System
TC28-1437 1.8ml
TC28-1438 2.2ml



Projection System
TC28-1439 1.8ml
TC28-1440 2.2ml

Шприцы стоматологические карпульные
Seringues à usage dentaire
Dentalspritzen

Шприці стоматологічні карпульні
Jeringas dentales
المحاقن الأسنان




Aluminium - Dark
TC28-1441 1.8 ml
TC28-1442 2.2 ml



Aluminium - Silver
TC28-1443 1.8 ml
TC28-1444 2.2 ml



Aluminium - Orange
TC28-1445 1.8 ml
TC28-1446 2.2 ml



Aluminium - Green
TC28-1447 1.8 ml
TC28-1448 2.2 ml



Aluminium - Blue
TC28-1449 1.8 ml
TC28-1450 2.2 ml

Костные кусачки
Hohlmeißelzangen
Pinces-gouges

Кусачки для кісток
Pinzas-gubias
رونجيورس العظام



TC28-1451
Blumenthal
15cm -90°



TC28-1452
Blumenthal
15cm -45°



TC28-1453
Blumenthal
15cm -30°



TC28-1454
Friedman
14cm



TC28-1455
Friedman 14.5cm
Thin
Fine



TC28-1456
Micro-Friedman 15cm
A form "S"
"S" Shaped



TC28-1457
Micro-Friedman
16cm



TC28-1458
Friedman
14 cm
30°



TC28-1459
Friedman
14 cm
45°



TC28-1460
Micro-Friedman
14cm -90°

Костные кусачки
Hohlmeißelzangen
Pincès-gouges

Кусачки для кісток
Pinzas-gubias
رونجيورس العظام



TC28-1461
Kleinert-Kulz
16cm



TC28-1462
Cohen 10cm
Tronchesini per gengive
Tissue nippers,
Pincès à couper le tissu
Alicatess para encías



TC28-1463
Cohen 14cm
Tronchesini per gengive
Tissue nippers,
Pincès à couper le tissu
Alicatess para encías



TC28-1464
Beyer
18.5cm




TC28-1465
Clevelanda
14cm

Костные кусачки
Hohlmeißelzangen
Pincès-gouges

Кусачки для кісток
Pinzas-gubias
رونجيورس العظام




TC28-1466
Rongeurs
16.5cm




TC28-1467
Rongeurs
Side Cutting
16.5cm



TC28-1468
Round beaks
Double Action
18 cm



TC28-1469
Ruskin-Mini
15.5 cm
Bone rongeurs,
multiple action



TC28-1470
Bohler
16 cm
Straight

Костные кусачки Кусачки для кісток
 Hohlmeißelzangen Pinzas-gubias
 Pincès-gouges روليجورس العظام



TC28-1471
Beyer
18 cm




TC28-1472
Böhler
15 cm
Bone cutting forceps,
multiple action



TC28-1473
Fig. 50
Liston
14cm



TC28-1474
Fig. 51
Liston
19cm



Kerrison
15 cm

TC28-1475
2 mm

TC28-1476
3 mm

TC28-1477
5 mm




TC28-1478
Tomati
12cm

- | | |
|--|--|
|  Костные кусачки |  Кусачки для кісток |
|  Hohlmeißelzangen |  Pinzas-gubias |
|  Pince-gouges |  رونجورس العظم |



TC28-1479
Universal




TC28-1480
Mead 1 A
16.5 cm



TC28-1481
Mead 2 A
16.5 cm



TC28-1482
Luer
14 cm



TC28-1483
Luer
14 cm




TC28-1484
Luer 15.5cm




TC28-1485
Luer, 17.5cm

Костные кусачки
Hohlmeißelzangen
Pincès-gouges

Кусачки для кісток
Pinzas-gubias
رونجيورس العظام



TC28-1486
Adson



TC28-1487
Cleveland 3
14 cm




TC28-1488
Cleveland 4
14 cm



TC28-1489
Cleveland 4 A
17 cm



TC28-1490
Cleveland 5
Side cutting
17 cm




TC28-1491
Cleveland 5 A
Side Cutting
17 cm



TC28-1492
Cleveland 5 S
15 cm
Side Cutting

Костные кусачки
Hohlmeißelzangen
Pincers-gouges

Кусачки для кісток
Pinzas-gubias
رونجيورس العظام




TC28-1493
Cleveland 10
17 cm




TC28-1494
Cleveland-Gardner
end Cutting



TC28-1495
Jansen
17.5 cm



TC28-1496
Jansen
18 cm



TC28-1497
Bane
Slight Curved
17.5cm



TC28-1498
Lempert Delicate
16cm

Инструменты для удаления коронок

Arrache-couronnes

Teleskopkronenzangen

Инструменти для видалення коронок

Arranca-corona

مزيلات التاج



TC28-1499
16cm

Crown Plier with plastic inserts Straight
For ceramical and plastic crowns and bridges, with
plastic inserts (sterilization up to 180° C).



TC28-1500
14.5cm

Crown Plier
with plastic inserts
Angled



TC28-1501
14cm

Universal plier for the removal
of ceramic crowns and bridges



TC28-1502
14cm

Universal plier for the removal
of ceramic crowns and bridges



TC28-1503
Furrer
Copper Ring
Removing Plier.
15cm
Curved



TC28-1504 12.5 cm
TC28-1505 15.0 cm
Böhm
Pliers for telescope
crowns.
Diamond points
coated.



TC28-1506
Nerve-Canal Pliers
straight



TC28-1507
Nerve-Canal Pliers
Curved



TC28-1508
Metal tip
22cm
Hammer With Wooden Handle



TC28-1509
Metal tip
22cm



TC28-1510
1
For plaster



TC28-1511
2
For metal



TC28-1512
Fret saw
25cm



TC28-1513
Small
Saw Blades




TC28-1514
Large
Saw Blades




TC28-1515
Ceri Saw



TC28-1516
Plaster Shears
20cm



Crucible Tongs
TC28-1517 25cm
TC28-1518 32cm



TC28-1519
Ceri Saw
18cm



TC28-1520
Blade
0,12x5 mm
Conf. 12 pcs



TC28-1521
Blade
0,25x5 mm
Conf. 12 pcs



TC28-1522
Blade
0,25x3 mm
Conf. 12 pcs



Wazirbad Road, Ugoki, Sialkot-51310 PAKISTAN

Phone: +92-52-3561347-3553588

Phone: +92-52-3561202

Cell: +92-300-6139615

Cell: +92-322-6139615

Email: info@tcmedicline.com

Email: exporttcmedicline@gmail.com